

# A-dec Wireless Module

| Product               | Model      |
|-----------------------|------------|
| A-dec Wireless Module | 43.0536.00 |

The A-dec Wireless Module device enables connection of the dental equipment to the A-dec Web App via Wi-Fi (Dual Band 802.11a/b/g/n/ac) and Bluetooth 5.0 BR/EDR/LE.



**NOTE** To satisfy RF exposure requirements, this device and its antenna must operate with a separation distance of at least 20 cm from all persons.

## Canada - ISED

| Product               | IC Number        |
|-----------------------|------------------|
| A-dec Wireless Module | 27025-ADEC430536 |

For indoor use only. This device complies with Industry Canada license-exempt RSS standard(s). Operation is subject to the following two conditions:

1. This device may not cause interference.
2. This device must accept any interference, including interference that may cause undesired operation of the device.

Utilisation en intérieur uniquement. Le présent appareil est conforme aux CNR d’ISDE Canada applicables aux appareils radio exempts de licence. L’exploitation est autorisée aux deux conditions suivantes :

1. l’appareil ne doit pas produire de brouillage, et
2. l’appareil doit accepter tout brouillage radioélectrique subi, même si le brouillage est susceptible d’en compromettre le fonctionnement.

## FCC ID

| Product               | FCC ID           |
|-----------------------|------------------|
| A-dec Wireless Module | 2AY33-ADEC430536 |

## FCC Compliance

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions:

1. This device may not cause harmful interference.
2. This device must accept any interference received, including interference that may cause undesired operation.

## FCC CAUTION

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

## Compliance with FCC Requirement 15.407(c)

Data transmission is always initiated by software, which is passed down through the MAC, through the digital and analog baseband, and finally to the RF chip. Several special packets are initiated by the MAC. These are the only ways the digital baseband portion will turn on the RF transmitter, which it then turns off at the end of the packet. Therefore, the transmitter will be on only while one of the aforementioned packets is being transmitted. In other words, this device automatically discontinues transmission in case of either absence of information to transmit or operational failure.

## Frequency Tolerance: $\pm 20$ ppm

This transmitter in this device must not be co-located or operated in conjunction with any other antenna or transmitter.



86.0924.00 Rev 1

Date of Issue: 2021-10-7

Copyright 2021 A-dec Inc.

All rights reserved.

IFUhalf



## **A-dec Headquarters**

2601 Crestview Drive  
Newberg, Oregon 97132  
United States  
Tel: 1.800.547.1883 within USA/CAN  
Tel: +1.503.538.7478 outside USA/CAN  
Fax: 1.503.538.0276  
[www.a-dec.com](http://www.a-dec.com)

## **A-dec Australia**

Unit 8  
5-9 Ricketty Street  
Mascot, NSW 2020  
Australia  
Tel: 1.800.225.010 within AUS  
Tel: +61.(0).2.8332.4000 outside AUS

## **A-dec China**

A-dec (Hangzhou) Dental Equipment Co., Ltd.  
528 Shunfeng Road  
Qianjiang Economic Development Zone  
Hangzhou 311100, Zhejiang, China  
Tel: 400.600.5434 within China  
Tel: +86.571.89026088 outside China

## **A-dec United Kingdom**

Austin House  
11 Liberty Way  
Nuneaton, Warwickshire CV11 6RZ  
England  
Tel: 0800.ADEC.UK (2332.85) within UK  
Tel: +44.(0).24.7635.0901 outside UK