

Label and Location

		Medical Machinery Product Name: Wrist blood pressure monitor Registrar Name: FITOLOGY Medical Technology (Hunan) Co., Ltd Registrar's Address: Room 017, Building 3, Medical and Health Industry Incubation Park, Nanxian Economic Development Zone, Nanchou Town, Nanxian County, Yiyang City, Hunan Province The Entrusted production enterprise name: Hunan YinuoSheng Medical Equipment Co.,Ltd Production Address: No.301,3rd Floor,Building 3,Emerging Industrial Park,Nanxian Economic Development Zone,Yiyang City,Hunan Province Input: 5V=1A Working Conditions: Temperature: -5℃~+45℃ Humidity: 10%RH-95%RH The condensation is not included; Atmospheric pressure: 80KPa-105kPa * This product has a built-in rechargeable lithium battery, battery use safety warning: disassembly, impact, extrusion or put into the fire is prohibited; If there is serious expansion, do not continue to use; do not put in a high-temperature environment; battery immersion is prohibited after the use of water. * Accessibility of some features may vary across regions. Product pictures, videos and display contents on the foregoing pages are provided for reference only. Actual product features and specifications (including but not limited to appearance, colour, and size), as well as actual display contents (including but not limited to backgrounds, UI, icons and videos) may vary. Consumer service hotline: 86-0571-86526367 Warning: 1.Do not wear the product during medical procedures, such as during magnetic resonance imaging (MRI) or the use of high-frequency medical devices. Keep the product away from any strong magnetic objects, or the product may fail. 2.The heart rate module feature on this product does not support asystole identification or ST segment detection. 3. In the event of arrhythmia (such as premature atrial or ventricular beats, or atrial fibrillation), bear in mind that the results may be inaccurate or inconclusive. 4. Diabetes, high cholesterol, and high blood pressure can affect peripheral circulation, and hence affect measurement results. 5. Blood pressure readings for users who suffer from blood circulation disorders or blood diseases are prone to error. 6. Avoid making newborns, infants who weigh less than 10 kg, or individuals who are unable to express themselves, wear or use this product.	

Label

Due to the product is too small to place such a large label, this label is placed on the package.

