

4. Label and Location

4.1 Label

[STERLINK U510]



Product name / Nom du produit : MD1-0194-02
Low temperature plasma sterilizer/ Stérilisateur plasma basse température
Brand name / Nom de marque : STERLINK™
Model name / Nom du modèle : STERLINK U510
Rating / Notation: 100-120VAC~, 50/60 Hz, 1.7kVA max.
Weight / Poids: 30 kg
FCC ID: 2BEXA-STU






Plasmapp Co.,Ltd.
A-1005, 10F, 670, Daewangpangyo-ro,
Bundang-gu, Seongnam-si, Gyeonggi-do,
13494, Republic of Korea

Please read the user manual carefully before use.
Veuillez lire attentivement le manuel d'utilisation avant utilisation.

This device complies with part 15 of the FCC Rules.
Operation is subject to the following two conditions:
(1) This device may not cause harmful interference, and
(2) this device must accept any interference received,
including interference that may cause undesired
operation.

Operation

Storage





www.plasmapp.com/eng/sub02-03.php

UDI



SN **STUA4ZD0001**  **YYYY-MM-DD**

(01)08800066500306
(11)250720
(21)STUA4ZD0001

[510]



Product name / Nom du produit : MD1-0287-00
Low temperature plasma sterilizer/ Stérilisateur plasma basse température
Brand name / Nom de marque : Plasma Steri
Model name / Nom du modèle : 510
Rating / Notation: 100-120VAC~, 50/60 Hz, 1.7kVA max.
Weight / Poids: 30 kg
FCC ID: 2BEXA-STU






Plasmapp Co.,Ltd.
A-1005, 10F, 670, Daewangpangyo-ro,
Bundang-gu, Seongnam-si, Gyeonggi-do,
13494, Republic of Korea

Please read the user manual carefully before use.
Veuillez lire attentivement le manuel d'utilisation avant utilisation.

This device complies with part 15 of the FCC Rules.
Operation is subject to the following two conditions:
(1) This device may not cause harmful interference, and
(2) this device must accept any interference received,
including interference that may cause undesired
operation.

Operation

Storage





www.plasmapp.com/eng/sub02-03.php

UDI



SN **STUA4ZD0001**  **YYYY-MM-DD**

(01)08800066500306
(11)250720
(21)STUA4ZD0001

4.2 Label Location

(단위: mm)



This label is permanently affixed to the back of the product and cannot be easily removed.