

LABEL AND LOCATION

Mode Overview

- Red & NIR: Comprehensive Repair Mode
- NIR: Deep Relief Mode (Invisible)
- Red: Skin Rejuvenation Mode
- Pulse: Energy Recovery Mode

Mode: Red + NIR / Pulse			Mode: NIR			Mode: Red		
6 inch (mW/cm ²)	12 inch (mW/cm ²)		6 inch (mW/cm ²)	12 inch (mW/cm ²)		6 inch (mW/cm ²)	12 inch (mW/cm ²)	
10W	23.3	17.1	5W	13.5	10.1	5W	6.6	4.6
25W	58	42.1	10W	29.3	21.8	10W	14.9	10.3
40W	94	68.3	20W	77.2	57.5	20W	32.3	22.2
55W	143.6	101.4	25W	90.6	67.6	30W	51.7	35.7

Note:

- Dose (J/cm²) = Irradiance (mW/cm²) × Time (s) × 0.001
- The actual dose depends on both the usage time and the irradiance from the product.
- 850nm Near-infrared(NIR) Light- Invisible but Active
- Pulse Mode: Constant 660nm Red Light + Pulsed 850nm NIR Light at 40Hz

Brand: Comfytemp
Model: K4037A
Input: 24V ~ 3A
Shenzhen Yicai Health Technology Co.,Ltd
Made in China
S/N:

This device complies with part 15 of the FCC Rules.
Operation is subject to the following two conditions:
(1) This device may not cause harmful interference, and
(2) This device must accept any interference received, including interference

FCC ID: 2BGBM-RLP

