





Description: Device label permanently affixed to the equipment	Material	Polyester
	Thickness	0,25 mm
	Code	L0-MCU-AB-00
	Date (mm/dd/yyyy)	03/12/2026

Device label position on the MostCare^{Up} Device

~ 100-240VAC
50/60Hz, 60VA
T1250mAh250V



REF
SN
FR:
UDI

MostCare^{Up}
Hemodynamic monitor

Bio-SiTM International srl
Via di Barbacane, 19
50133 (FI) Italy
www.bio-si.it

FCC ID: 2BOMM-00C01
Contains FCC ID: QQWF111

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

MD
Rx Only
IP22
ITA

Device label position on the MostCare^{Up} Device

