

Xlink 2.0

Wireless Monitor Bridge

User Manual


Scan or visit to get started:



xchart.com/xlink

Dongle


Front Side



- ① FIND PAIR Button: Press to identify the paired adapter.
- ② USB-C Connector: Used to plug the dongle into your computer.
- ③ ACTIVITY LED: Indicates when the paired adapter is within range.


Adapter

Front Side



- ① CONFIG Button: Press to enter configuration mode.
- ② POWER Connector: Supplies power to the adapter.
- ③ ETHERNET Connector: Used to connect the adapter to a medical device via Ethernet.

Back Side



- ④ PAIRED DONGLE PRESSED: Glows if FIND PAIR button pressed.
- ⑤ ETHERNET MODE: Glows when configured for Ethernet mode.
- ⑥ SERIAL MODE: Glows when configured for Serial mode.
- ⑦ DATA ACTIVITY: Blinks when sending or receiving data.
- ⑧ SERIAL Connector: Used to connects the adapter to a medical device via a Serial port.

FCC Compliance Statement

This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions:

This device may not cause harmful interference, and

This device must accept any interference received, including interference that may cause undesired operation.

This equipment has been tested and found to comply with the limits for a Class B digital device and an intentional radiator, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation.

This equipment generates, uses, and can radiate radio frequency energy, and if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation.

If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

FCC RF Exposure Information and Statement The SAR limit of USA (FCC) is 1.6 W/kg averaged over one gram of tissue.

Device types: It has also been tested against this SAR limit. This device was tested for typical body-worn operations with the back of the product kept 5mm from the body. To maintain compliance with FCC RF exposure requirements, use accessories that maintain a 5mm separation distance between the user's body and the back of the product. The use of belt clips, holsters and similar accessories should not contain metallic components in its assembly. The use of accessories that do not satisfy these requirements may not comply with FCC RF exposure requirements, and should be avoided. Body-worn

Operation This device was tested for typical body-worn operations. To comply with RF exposure requirements, a minimum separation distance of 5mm must be maintained between the user's body and the product, including the antenna. Third-party belt-clips, holsters, and similar accessories used by this device should not contain any metallic components. Body-worn accessories that do not meet these requirements may not comply with RF exposure requirements and should be avoided. Use only the supplied or an approved antenna.