


Diasorin Liaison XS RFID System Description

2025-11-21

12. External Photographs

1. Reagent Loading Bay



- The four RFID antennas for the max. 4 inserted ancillary packs are located at the left side behind a plastic wall to read/write information on the RFID tag attached to the ancillary pack.
- The ten RFID antennas for the max. 10 inserted reagent racks are located below the 10 lanes to read/write information on the RFID tag attached to the reagent racks.

2. Starter Reagent Module



- The two RFID antennas for the starter reagent racks are located behind the rear panel to read/write information on the RFID tags attached to the starter reagent bottles.