

Differences declaration

To whom concern:

Our company Wuhan Youkey Bio-Medical Electronics Co., Ltd.

Here by declares:


Multi-model differences are as follows:

Main model: Listed model:		Main Model No.: Q70				
		Schematic diagram 电路原理图	PCB layout PCB 走线	Appearance 外观结构	Model name 型号命名	Other 其它
Series	Model No.:	☒	☒	△	△	☒
Q31,Q32,Q33,Q34,Q35,Q36,Q37,Q38,Q39,Q39L (Sono-Key),Q80,Q51,Q52,Q53,D80,Q7,MyLabPocket						

Note:

- If it is a same, mark it as ☒
如相同, 则标识为 ☒
- If there is a difference, mark it as △; If there are differences other than model name differences, please describe them in detail
如有差异, 则标注为△; 除了型号命名差异外, 如有其它差异则需要下面详细描述

Differences describe
Our company's research, development, and production of the Pocket Ultrasound System cover models Q7, Q70, MyLab Pocket, Q31, Q32, Q33, Q34, Q35, Q36, Q37, Q38, Q39,Q39L (Sono-Key),D80, ,Q80,Q51, Q52, and Q53. The devices of these models utilize identical hardware circuitry, transducers, software, and algorithms. Q70 serves as the primary tested model. Compared to the Q70, the Q7 and MyLab Pocket feature the same types of transducers, which are interchangeable, and differ only in their external appearance and labeling. In contrast, when compared to the Q7, the models Q31, Q32, Q33, Q34, Q35, Q36, Q37, Q38, Q39,Q39L (Sono-Key),D80, ,Q80,Q51, Q52, and Q53 are equipped with non-interchangeable transducers and also feature different appearances and labeling. The only difference on the label is the model.

Name:	<u>zeyi hu</u>
Date:	<u>2026-07-17</u>
Signature:	<u></u>
Title:	<u>Corporate juridical person</u>
Phone number:	<u>+86 18971511819</u>
Email:	<u>huzeyi@youkeymedical.com</u>