



Radio Regulatory Compliance Information

MyCareLink Home™ Communicator Model 26950

FCC compliance information

FCC ID: LF526950

Contains FCC ID: LF526900

This device complies with Part 15 of the FCC rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

The user is cautioned that changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

For radiofrequency exposure safety, the minimum separation distance between the MyCareLink Home base and the patient's body is 20 cm. The cabled telemetry reader located inside the MyCareLink Home base can be used without minimum separation distance limitations.



Medtronic, Inc.

710 Medtronic Parkway
Minneapolis, MN 55432
USA

www.medtronic.com

+1 763 514 4000



M052868C001

© 2025 Medtronic
M052868C001 A
2025-04-02