

FCC E-Label Information

302-000-000 McGRATH MAC + Connected Video Laryngoscope

1. Users can access the required information on the electronic screen of the McGRATH MAC + Connected Video Laryngoscope in no more than three steps from the main menu. The steps are:

To access FCC ID follow these steps:

Step 1) With the laryngoscope powered off, press and hold the power button to access the "Settings" Menu.

Step 2) Continuing to hold the power button, select "User Manual".

No special access codes or permissions are required to go through the above steps to access the FCC and IC IDs.

2. Image of the e-label screen is provided below:

