


## Care Technologies AC3 Transmitter Users Guide



The AlertCare3 transmitter is for use in conjunction with the Care Technologies Community Care System. To use the transmitter, simply affix it to the resident using the supplied permanent snap band. The transmitter can be affixed to the resident's wrist or ankle.

It is possible that residents in wheelchairs may not be able to be monitored using the AC3 transmitter. The metal frame of the wheelchair may block the signal sent out by the transmitter.

FCC NOTICE for AC3 Transmitter FCC ID: QCKAC3TX

This device complies with Part 15 of the FCC rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Changes or modifications not specifically authorized by Care Technologies could void the user's authority to operate the equipment and void any warranty expressed or implied.