

REV	ECO #	DATE	DESCRIPTION	DWG	APRV
XX	0XXX	DD MMM YYYY	DES OF CHG	XXX	

American TeleCare

NX Patient Station


www.americantelecare.com

NX Patient Station FCC ID: QQE206A

Voltage: 12 V

Modem FCC ID: LNQUSA-31455-MME

Current: 2.4 A Max



DANGER: RISK OF EXPLOSION

IF USED IN THE PRESENCE OF


FLAMMABLE ANESTHETICS.

Grounding reliability can only be achieved

when this equipment is connected to

a receptacle marked "Hospital Grade".

NX Patient Station




Tested to Comply

With FCC Standard

FOR HOME OR OFFICE USE

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.



This document contains confidential and proprietary information of American TeleCare, Inc. and cannot be used, duplicated, or disclosed without the express written permission of American TeleCare, Inc.



American TeleCare, Inc.

7640 Golden Triangle Drive

Eden Prairie, MN 55344

952-897-0000


www.americantelecare.com

TOLERANCES

FRAC	ANG	DEC	TOL
1/X	X°	.X	±.1
1/XX	X.X°	.XX	±.01
	X.XX°	.XXX	±.005
		.XXXX	±.0005

ANSI Y14.5-1994

THIRD ANGLE PROJECTION



ENGINEER

JOHN BLOMBERG

DD MON YYYY

DRAFTER

JOHN BLOMBERG

DD MON YYYY

APPROVED

C500-40110-MD

REV C

SIZE D

SCALE 8:1

SHEET 1 OF 1

File Name: C500-40110 REV C p2 Date: Friday, January 31, 2003