

American TeleCare, Inc.
Model: NX Patient Station 206

02077-02	<u>Wi</u>	<u>P206-003-00</u>	<u>B</u>	<u>6/18/02</u>			
PartNumber	Qty.				RefDes		
PW-7391-0001	1	BD PUR PW-7391-0001 REV A	PUR AS SINGLE BOARD		PWB		
						[EQV]	
SR-0301-4100	27	SMR 0805 10.0K 1% 1/10W	100PPM	805	R201*,202*,303*,309*,311*,312*,506*,507*,513*,810*,		
						[EQV]	DALE PANASONIC CRCW08051002F ERJ6ENF1002V
SR-0305-0000	8	SMR 0805 0 5% 1/10W	200PPM	805	R308*,313*,512*,514*,606*		
						[EQV]	ROHM KOA MCR10EZHMJW000 RM73Z2A
SR-0301-3100	2	SMR 0805 1.00K 1% 1/10W	100PPM	805	R310*,901*		
						[EQV]	PANASONIC ROHM ERJ6ENF1001V MCR10 EZHM FX 1001
SR-0305-1220	2	SMR 0805 22 5% 1/10W	200PPM	805	R314*,515*		
						[EQV]	DALE ROHM CRCW0805220J MCR10EZHMJW220
SR-0301-5100	20	SMR 0805 100K 1% 1/10W	100PPM	805	R401*,402*,403*,404*,405*		
						[EQV]	KOA PANASONIC RK73H2A1003F ERJ6ENF1003V
SR-0305-5220	1	SMR 0805 220K 5% 1/10W	200PPM	805	R601*		
						[EQV]	KAMAYA KOA RMC10-224JT RM73B2A224J
SR-0336-0002	2	SMR 2010 0.02 OHM 5% 1W	100 PPM	2010	R602*,603		
						[EQV]	IRC IRC LRC-LRF2010-01-R020-J LRC-LRF2010-01-R020-
SR-0301-5274	1	SMR 0805 274K 1% 1/10W	100PPM	805	R604*		
						[EQV]	ROHM PANASONIC MCR10EZHMF2743 ERJ6ENF2743V
SR-0301-5115	1	SMR 0805 115K 1% 1/10W	100PPM	805	R605*		
						[EQV]	KAMAYA RMC10-1153F
SR-0305-1100	1	SMR 0805 10 5% 1/10W	200PPM	805	R608*		
						[EQV]	KAMAYA ROHM RMC10-10RJT MCR10EZHJW100
SR-0301-4750	1	SMR 0805 75K 1% 1/10W	100PPM	805	R716*		
						[EQV]	KAMAYA RMC10-7502F
SR-0305-3300	1	SMR 0805 3.0K 5% 1/10W	200PPM	805	R717*		

American TeleCare, Inc.
Model: NX Patient Station 206

						[EQV]	KAMAYA	RMC10-302JT
SR-0305-4150	8	SMR 0805 15K 5% 1/10W	200PPM	805	R803*,804*,805*,806*,807*			
						[EQV]	KAMAYA	RMC10-153JT
SR-0305-3150	3	SMR 0805 1.5K 5% 1/10W	200PPM	805	R809*,833*,903*			
						[EQV]	ROHM	MCR10EZHMJW152
							KOA	RM73B2A152J
SR-0305-5180	1	SMR 0805 180K 5% 1/10W	200PPM	805	R902*			
						[EQV]	KAMAYA	RMC10-184J
SC-0421-6220	4	1206 X7R 16V 10% 2.2UF	SM CAP CERAMIC	1206	C201*,202*,203*,915*			
						[EQV]	CAL CHIP	GMC31X7R225K16NT
							PANASONIC	ECJ-3YB1C225K
SC-0405-3100	37	0805 X7R 50V 10% 1000PF	SM CERAMIC CAP	805	C201A*,201B*,201C*,201D*,201E*,201F*,201G*,201H			
						[EQV]	NIC	NMC0805X7R102K050
SC-0405-5100	22	0805 X7R 50V 10% .1UF	SM CERAMIC CAP	805	C204*,501*,502*,503*,504*			
						[EQV]	VITRAMON	VJ0805Y104KXAMT
							KEMET	C0805C104K5RAC
SC-0500-2100	4	0805 COG 50V 5% 100PF	SM CERAMIC CAP	805	C510*,511*,512*,513*			
						[EQV]	MURATA ERIE	GRM40COG101J050
							MURATA	GRM40COG101J100AD
SC-0405-4100	5	0805 X7R 50V 10% .01UF	SM CERAMIC CAP	805	C601*,908*,909*,911*,913*			
						[EQV]	MURATA	GRM40X7R103K050A
							MURATA	GRM40X7R103J050AL
SC-0409-6100	2	0805 X7R 16V 10% 1UF	SM CAPACITOR CERAMIC	805	C602*,626*			
						[EQV]	CAL CHIP	GMC21X7R105K16NT
							KEMET	C0805C105K4RAC
SC-0243-8221	2	7343 TANT 10V 10% 220UF	SM CAP MOLDED EXTD RANGE	7343	C605*,606*			
						[EQV]	KEMET	T491D227K010AS
							NEMCO	PCT220/10DK
SC-0214-7100	10	6032 TANT 16V 10% 10UF	SM CAPACITOR MOLDED	6032	C607*,620*,621*,622*,623*			
						[EQV]	PHILIPS	49MC106C016KOASFT
							KEMET	T491C106K016AS
SC-0404-5220	1	0805 X7R 25V 10% .22UF	SM CERAMIC CAP	805	C608*			
						[EQV]	KEMET	C0805C224K3RAC
							AVX	08053C224KAT
SC-0500-1470	4	0805 COG 50V 5% 47PF	SM CERAMIC CAP	805	C610*,903*,910*,912*			
						[EQV]	PANASONIC	ECUV1H470JCG
							KEMET	C0805C470J1GAC

American TeleCare, Inc.
Model: NX Patient Station 206

SC-5375-7470	1	7343 TANT 16V 20% 47UF	SM CAP MOLDED LOW ESR	7343	C611*			
						[EQV]	AVX	TPSD476M016R0150
SC-4042-7100	1	ECE-V1VA100SR PANA C-SIZE	SM CAP AL EL 10UF 20% 35V	SMD	C616			
						[EQV]	PANASONIC	ECE-V1VA100SR
SC-5327-6470	2	6032 TANT LOW ESR 4.7UF	SM CAP MOLDED 35V 20%	6032	C618*,619*			
						[EQV]	AVX	TPSC475M035R0600
SC-0500-1270	2	0805 COG 50V 5% 27PF	SM CERAMIC CAP	805	C801*,802*			
						[EQV]	MURATA ERIE KEMET	GRM40COG270J050B C0805C270J5GAC
SC-4040-7330	4	ECE-V1VA330WP PANASONIC	SM CAP AL EL 33UF 20% 35V	SMD	C843,844*,845*,846			
						[EQV]	PANASONIC	ECE-V1VA330WP
SC-1100-0180	2	0805 COG 50V 1.8PF+-.25PF	SM CERAMIC CAP	805	C901*,902*			
						[EQV]	MURATA ERIE MURATA ERIE	GRM40COG1R8C100A GRM40COG1R8C050B
SC-0500-2220	1	0805 COG 50V 5% 220PF	SM CERAMIC CAP	805	C904*			
						[EQV]	MURATA	GRM40COG221J050BD
SC-1200-1100	1	0805 COG 10PF +-.5PF 50V	SM CERAMIC CAP	805	C906*			
						[EQV]	AVX MURATA ERIE	08055A100JAT2A GRM40COG100C050BL
SC-1100-1100	1	0805 COG 10PF.+-.25PF 50V	SM CERAMIC CAP	805	ALT C906*			
						[EQV]	MURATA	GRM40COG100C050
SC-1200-0330	1	0805 COG 3.3PF +-.5PF 50V	SM CERAMIC CAP	805	C907*			
						[EQV]	VITRAMON	VJ0805A3R3DXA
SC-1100-0330	1	0805 COG 3.3PF+-.25PF 50V	SM CERAMIC CAP	805	ALT C907*			
						[EQV]	MURATA ERIE	GRM40COG3R3C050BD
SD-1709-1200	2	MURS120T3 MOTO SMB	DIODE RECTIFIER 1AMP 200V	SMB	D601*,602*			
							MOTOROLA	MURS120T3
SD-1000-0001	5	PMBD 914/4148 SOT23	SM SINGLE DIODE GEN PURP	SOT23	D603*,604*,605*,606*,607*			
						[EQV]	ALLEGRO MOTOROLA	TMPD914LT MMBD914LT1
ST-8709-4030	4	IRF7403 INTERNATIONAL	HEXFET POWER MOSFET SON8	SON8-150	Q601,602,603,604			

American TeleCare, Inc.
Model: NX Patient Station 206

						[EQV]	INTERNATIONAL R	IRF7403
SS-5000-1102-01	1	SA1110BD M/B:INTEL 256BGA	MICROPROCESSOR STRONGARM	256BGA	U201			
						[M/B]	INTEL	SA1110BD
							INTEL	GDS1110BD
SS-5000-1101-01	1	GDS1111-BA M/B:INT 256BGA	MICROPROCESSOR STRONGARM	256BGA	U202			
						[M/B]	INTEL	GDS1111-BA
SS-1401-7450-01	1	SN74LVC32245GKER M/B:TI	TRANSCEIVER 32BIT 96BGA	96BGA	U203			
						[M/B]	TI	SN74LVC32245GKER
SS-3100-0280-01	2	RC28F128J3A-150 M/B:INTEL	STRAFLASH 64BGA UNPROGRAM	64BGA	FOR U301,302 UNPROGRAMMED			
						[M/B]	INTEL	RC28F128J3A-150
SS-3100-0280-50	2	PROGRAM SS-3100-0280-01	STRAFLASH MEMORY 64BGA	64BGA	U301,302 PROGRAMMED			
						[M/B]	ORDER PER PRINT	SS-3100-0280-50
SS-3200-6320-01	2	K4S561632C-TC7C M/B:SAMU	SDRAM CMOS 54TSOP 133MHZ	54TSOP	U303*,304*			
						[M/B]	SAMSUNG	K4S561632C-TC7C
SS-2502-1600-01	1	MAX1600EAI M/B:MAXIM	POWER SWITCHING 28SSOP	SSOP-28	U404			
						[M/B]	MAXIM	MAX1600EAI
SS-1401-3340-01	1	LTC1334CG M/B:LIN 28SSOP	TRANSCEIVER 5V RS232/485	SSOP-28	U501			
						[M/B]	LINEAR TECH	LTC1334CG
SS-1061-3050-01	1	DS1305E M/B:DALLAS	REAL-TIME CLOCK 20TSSOP	TSSOP-20	U504			
						[M/B]	DALLAS	DS1305E
							DALLAS	DS1305E
SS-1264-1300-01	1	UCB1300BE M/B:PHILIPS	AUDIO/TELECOM CODEC QFP48	QFP48	U505			
						[M/B]	PHILIPS	UCB1300BE
SS-1382-5750-01	1	MAX4575EUA M/B:MAX 8UMAX	SWITCH ANALOG, DUAL SPST	8uMAX	U507			
						[M/B]	MAXIM	MAX4575EUA
SS-2502-6310-01	1	MAX1631CAI M/B:MAX 28SSOP	POWER-SUPPLY CONTROLLER	SSOP-28	U601			
						[M/B]	MAXIM	MAX1631CAI
SS-2562-6920-01	1	MAX1692EUB M/B:MAX 10UMAX	PWM REGULATOR STEP-DOWN	10uMAX	U602			
						[M/B]	MAXIM	MAX1692EUB
SS-2001-4119-01	1	MC34119D M/B MOTO	AUDIO AMP. SON8-150	SON8-150	U702			
						[M/B]	MOTOROLA	MC34119D

American TeleCare, Inc.
Model: NX Patient Station 206

SS-9999-0042-01	1	TUSB2077APT M/B:TI 48QFP	HUB 7-PORT FOR USB	QFP48	U801			
						[M/B]	TI	TUSB2077APT
SS-2502-0440-01	1	TPS2044ADR M/B:TI SON16	PWR DISTRIBUTION SWITCH	SON16-150	U802			
						[M/B]	TI	TPS2044ADR
SS-5000-0606-01	2	SN75240PW M/B:TI TSSOP-8	TRANSIENT VOLT SUPPRESSOR	TSSOP-8	U804,805			
						[M/B]	TEXAS INST.	SN75240PW
SS-9999-0044-01	1	BC212015BDN-E4 M/B:CSR	BLUETOOTH SYSTEM BGA-96	BGA96	U901			
						[M/B]	CSR	BC212015BDN-E4
SS-3100-3980-01	1	SST39VF800-90-4C-B3K M/B	FLASH MEMORY 8M BGA48	BGA48	U902			
						[M/B]	SST SST	SST39VF800-90-4C-B3K SST39VF800A-90-4C-
SS-9999-0045-01	1	CAF95901 M/B:CENTURION	SM ANTENNA BLUECHIP	SMD	U903			
						[M/B]	CENTURION	CAF95901
SL-4605-2450-01	1	MDR741F M/B:SOSHIN	FILTER BANDPASS 2.45GHZ	SMD	U904			
						[M/B]	SOSHIN	MDR741F
SL-6500-0622-01	1	NL322522T-220J M/B:TDK	SM INDUCTOR 22UH 5% 1210	1210	L201*			
							TDK	NL1210-220J
SB-2254-2600-01	3	DA1206D600R-00 M/B:STEWAR	SM FERRITE CHIP ARRAY	SMD	LC501*,502*,503*			
						[M/B]	STEWARD	DA1206D600R-00
SL-4611-0610-01	2	DO3316P-103 M/B:COILCRAFT	SM INDUCTOR PWR 10UH 20%	SMD	L601,602			
						[M/B]	COILCRAFT	DO3316P-103
SL-4614-0610-01	1	DO1608C-103 M/B:COILCRAFT	SM INDUCTOR PWR 10UH 20%	SMD	L603			
						[M/B]	COILCRFAT	DO1608C-103
SL-4801-0800-01	1	CM5018R800R-00 M/B:STEWAR	SM CHOKE COMMON MODE SO16	SMD	LC801			
						[M/B]	STEWARD	CM5018R800R-00
SB-2254-2300-01	2	DA1206E300R-00 M/B:STEWAR	SM FERRITE CHIP ARRAY	SMD	LC802*,803			
						[M/B]	STEWARD	DA1206E300R-00
SL-9013-0239-01	4	LL2012-FH3N9K M/B:TOKO	SM INDUCTOR 3.9NH 0805	805	L901*,902*,903*,904*			
						[M/B]	TOKO TOKO	LL2012-F3N9K LL2012-FH3N9K
SL-9900-0016-01	1	HHM1517 M/B:TDK	SM TRANSFORMER CHIP	SMD	T901			
						[M/B]	TDK	HHM1517
SY-2318-0363-01	1	HC49SD-3.6864CBF1601 M/B	SM CRYSTAL 3.6864MHZ FOX	SMD	X201*			
						[M/B]	FOX	HC49SD-

American TeleCare, Inc.
Model: NX Patient Station 206

SY-2371-0364-01	1	FOXSD/0368-20 M/B:FOX	SM CRYSTAL 3.6864MHZ 20PF	SMD	ALT X201*			
						[EQV]	FOX	FOXSD/0368-20
SY-2341-0015-01	2	FSR327 M/B:FOX 32.768KHZ	SM CRYSTAL WATCH 12.5PF	SMD	X202*,501			
						[M/B]	FOX ELECTRONICS	FSR327
SY-2380-1600-01	1	TSX-10A-16MHZ M/B:TOYOCOM	SM CRYSTAL 16.0000MHZ	SMD	X901			
						[M/B]	TOYOCOM	TSX-10A-16MHZ
SH-2015-0018-01	1	04-6214-018-000-800 M/B:A	SM CONN ZIF VERT 18POS	SMT-18POS	J502			
						[M/B]	AVX	04-6214-018-000-800
SH-2016-0018-01	1	52559-1890 M/B:MOLEX	SM CONN ZIF VERT 18POS	SMT-18POS	ALT J502			
						[M/B]	MOLEX	52559-1890
SH-2333-0006-01	1	557314-1 M/B:AMP 6POS RA	SM MODULAR JACK W/PEGS	SMT-6POS	J504*			
						[M/B]	TYCO ELECTRONIC	557314-1
							AMP	557314-1
SH-9006-0010-01	1	TSM-105-01-T-DV-P M/B:SAM	SM TERM STRIP 2X5 TIN P&P	SMT-2X5	J505			
						[M/B]	SAMTEC	TSM-105-01-T-DV-P
SA-1001-0002	1	3000 KEYSTONE 12MM DIA	SM COIN CELL RETAINER	SMD	J601			
						[EQV]	KEYSTONE	3000
SH-9005-1020	1	TSM-110-01-L-DV-M SAMTEC	SM TERM STRIP 2X10 10U"AU	SMT-2X10	J901			
						[EQV]	SAMTEC	TSM-110-01-L-DV-M
SH-2701-0009-01	1	K202HT-E9P-N M/B:KYCON 9M	SM D-SUB CONN W/BD LCK RA	SMT-DSUB9	P503			
						[M/B]	KYCON	K202HT-E9P-N
PS-0440-0813	2	7228-5 KEYSTONE 4-40X1/4L	SCREW JACK DSUB STEEL ZNC	MISC	PS5,6			
						[EQV]	KEYSTONE	7228-5
SH-1512-0002-01	1	53398-0290 M/B:MOLEX 1X2	SM HEADER W/SHROUD TIN	SMT-1X2	P702			
						[M/B]	MOLEX	53398-0290
LP-1016-4200-01	1	3306K-1-203 M/B BOURNS	TRIMPOT 20K RA NO-WASH	LEADED	R508			
						[M/B]	BOURNS	3306K-1-203
LY-2377-0601-01	1	FOX060 M/B:FOX HC49U	CRYSTAL 6.000MHZ SERIES	LEADED	X802			
						[M/B]	FOX ELECTRONICS	FOX060
NC-1304-0001	1	CI-192-028 EQV:BIVAR	SPACER, CRYSTAL, MYLAR	MISC.	SP1			
						[EQV]	BIVAR	CI-192-028

American TeleCare, Inc.
Model: NX Patient Station 206

CC-2053-0068-01	2	IC1FB-68PD-1.27DS-EJ M/B	CONN PCMCIA 68POS HIROSE	68 POS	J401,402*			
						[M/B]	HIROSE	IC1FB-68PD-1.27DS-EJ
MS-0620-0809	8	PHMS 6-20X1/4 STL/ZINC	SCR TYP AB TAP PAN HD PHL	MISC	MS5,6,7,8,9*,10*,11*,12*			
						[EQV]	PATTON	13N25XP0Z/AB
CC-2236-0004-01	1	558872-1 M/B:AMP 4POS	MODULAR JACK W/PEGS	LEADED	J503*			
						[M/B]	AMP	558872-1
CC-8007-0004-01	1	71226-0421 M/B:MOLEX 4POS	CONNECTOR ZIF RECEPT VERT	4 POS	J506			
						[M/B]	MOLEX	71226-0421
CC-1500-0022-01	1	PC722A M/B:SWITCHCRAFT	POWER JACK DC 3 TERMINAL	LEADED	J602*			
						[M/B]	SWITCHCRAFT	PC722A
CC-2227-0004-01	2	KUSBV-AS-1-N-WHT M/B:KYCN	CONN USB JACK VERT 4POS	4 POS	J801*,802*			
						[M/B]	KYCON	KUSBV-AS-1-N-BLK
							KYCON	KUSBV-AS-1-N-WHT
CC-1483-0008-01	1	787617-1 M/B:AMP 8POS RA	RECEPTACLE W/BOARD LOCK	8 POS	J803			
						[M/B]	AMP	787617-1
TP-9354-0008-01	1	SSW-104-02-T-Q-RA M/B:SAM	SOCKET STRIP RA 2X4 TIN	2 X 4	P504			
						[M/B]	SAMTEC	SSW-104-02-T-Q-RA
LS-2064-6160-01	1	8M122616 M/B:ERG	INVERTER DC-AC SUB ASSY.		SUB1			
						[M/B]	ERG	8M122616
PS-9002-4187	2	SPACER NYLON	.25 OD X .12 ID X .187 LN		PS1,2			
						[EQV]	RICHCO	R911-3
							BIVAR	9911-187
MS-0440-1231	2	PHMS 4-40X3/8" STAIN STL	SCREW PAN HEAD PHILIPS		MS1,2			
						[EQV]	ADAMS	94400
							PATTON	MS-0440-1231
MN-0440-0001	2	NUT HEX 4-40 UNC MACH SCR	STAINLESS STEEL		MN1*,2*			
						[EQV]	PATTON	4-40 HEX M/S NUT SS
							ADAMS	4-40 HEX M/S NUT SS
							ADAMS NUT& BOLT	87000
OMIT	48				R509*,510*,511*,701*,703*,704*,705*,706*,707*,708*,			
						[EQV]		
OMIT					U503,701; X701;			
						[EQV]		

American TeleCare, Inc.
Model: NX Patient Station 206

					[EQV]		
--	--	--	--	--	-------	--	--