





NOTES:

- 1. JOIN (3) TO (2) BY SONIC WELDING.
- 2. ASSEMBLE (1) TO (2) & (3) AFTER STEP 1 IS COMPLETE.
- 3. REFER TO PROE MODELS: 3983-000.ASM, 3983-MST.PRT.

REVISION HISTORY				
REV	EN	DESCRIPTION	DATE	EN ORIGINATOR
01	--	INITIAL ENG RELEASE	4-21-03	--
02	--	ADDED LABELS, ITEMS 5 - 9	6-12-03	--
03	--	REMOVED: 4194-000 & 4355-000. RENUMBERED ITEMS	9-10-03	--

UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN INCHES TOLERANCES ON: .X ± 0.1 .XX ± 0.01 .XXX ± 0.001 ANGLES ± 0.5	DRAWN BY ADB	DATE 4-21-03	 NONIN MEDICAL, INC. PLYMOUTH, MN The information disclosed herein is confidential and is proprietary to Nonin Medical, Inc. Nonin reserves all patent, proprietary and reproduction rights therein. This information is not to be used or disclosed to any persons without the express written authorization of a representative of Nonin Medical, Inc.		CONFIDENTIAL	
	CHECKED BY	DATE				
	ELECTRONIC FILE NAME 3983-000.DRW		TITLE FINAL ASSY, PATIENT OXIMETER, BLUETOOTH			
	DRAWING CONFORMS TO ANSI Y 14.5M-1994		SIZE B	DWG NO 3983-000	REV 03	
THIRD ANGLE PROJECTION 	DO NOT SCALE DRAWING		SCALE N/A	SHEET 1	OF 1	