

Label:



This device complies with Part 15 of the FCC Rules. Operation is subject to the following conditions: (1) This device may not cause harmful interference, (2) this device must accept any interference received, including interference that may cause undesired operation.

FCC ID: SMQ-LX9EDAN

Label Location

EDAN Diagnostic Ultrasound System  0123

Model: Acclarix LX9  2021-08-27

SN  560595-M21813020004 

Anti-electro-shock Type: Class I Max. Mass: 83kg

100V-240V~ 50Hz/60Hz 5.2A-2.5A

EC REP Shanghai International Holding Corp. GmbH

Eiffelstrasse 80, 20537 Hamburg, Germany

Tel: +49-40-2513175 E-mail: shholding@hotmail.com Rx Only

Edan Instruments, Inc.
#15 Jinhui Road, Jinsha Community,
Kongzi Sub-District, Pingshan District,
518122 Shenzhen, P.R. China

Tel: +86-755-26898326 (010)06844413811393
Fax: +86-755-26898330 (21)560595-M21813020004
http://www.edan.com (11)210827

VIDEO  DVI 

USB  VGA 

USB  RJ45  HDMI 

1

2

3

4

5

6

7

8

REV.

ECN

DESCRIPTION

	DRAW
--	------

CHECK

	APPR.
--	-------

DATE

1.1

SMOLX9EDAN更新为SMQ-LX9EDAN

[illegible]

1

[illegible]

Text is 5pt.

Color is Warm Grey 8

Full Round on each end

9mm


68mm

This device complies with Part 15 of the FCC Rules. Operation is subject to the following conditions: (1) This device may not cause harmful interference, (2) this device must accept any interference received, including interference that may cause undesired operation.

FCC ID: SMQ-LX9EDAN

技术要求

- 1.丝印质量和丝印检验规范请参考理邦公司《面贴、保护屏检验作业指导书》执行；
- 2.带公差标注的尺寸为重点检验和控制的尺寸；
- 3.未注公差按SJ/T 10628-1995-7级执行；
- 4.各图案所表示的功能请参考图面说明内容；
- 5.零件需符合《理邦有害物质管控标准》；
- 6:文档中相同色块区域表示相同加工工艺。

DESIGN	刘 期 盼	TITLE L X9 FCC声明标贴	PRODUCT Acclarix Lx9			
CHECK			DOC NO. EDAN\WI\6-01.55.469695	REV. 1.1	UNIT mm	SCALE 1:1
ENG.		MATERIAL 合成纸(哑膜)	SHEET 1 OF 1		PROJECTION 	
APPROVAL			 EDAN 深圳市理邦精密仪器股份有限公司 EDAN INSTRUMENTS, INC.			