





Check the image below for the correct details:


1. Select the Settings button on the main screen.



2. Select Device Information on the settings screen.




3. Check the FCC ID or IC number in the Device Information menu.




1.6.1 Locations of Labels

(1) XD4000 PX



“XD4000 PX” is listed on the rating label below.





Note1. Access to regulatory information is provided through the standard user interface and does not require any special codes, accessories, or permissions beyond the normal screen unlock or user sign-in process.


Note2. The FCC ID & IC Number is displayed electronically in a clear and legible manner and can be read without the aid of magnification

Note3. The FCC ID & IC Number is programmed by the responsible party and is secured such that it cannot be modified by third parties.

The regulatory information is implemented in a factory-set, unalterable format and is presented through a dedicated and exclusive application.

The FCC ID & IC Number is retrieved and displayed locally on the device and is accessible under the control of the end user in possession of the product.

The method is secure and ensures that the regulatory information is displayed correctly



XD4000 PORTABLE E

XD4000 PX

SN XDYJE-0001

UDI



(01)04547410564136
(11)251017
(21)XDYJE-0001



2025-10-17

FCC ID : W2Z-01000017
IC : 7736B-01000017



FUJIFILM Corporation

26-30, NISHIAZABU 2-CHOME, MINATO-KU, TOKYO 106-8620, JAPAN

Portable X-ray Unit FDR Xair II

EC REP

FUJIFILM Healthcare Europe GmbH
Balcke-Duerr-Allee 6, 40882 Ratingen, Germany



<https://med-lib.fujifilm.com/>

CH REP

FUJIFILM Healthcare Europe Holding AG
Sumpfstrasse 15, 6312 Steinhausen, Switzerland
Imported to Australia by: FUJIFILM Australia Pty Ltd
54 Waterloo Road, Macquarie Park, NSW 2113, Australia

CE 0197

MADE IN KOREA