


Regulatory compliance

To access regulatory compliance information on the INX 660 microscope, select Main Menu > System Settings > System Info > Regulatory Info.

Figure 1-2: Accessing regulatory compliance information



There are no special codes, accessories or permissions required to access the regulatory information with the exception of standard security protection to unlock the screen per KDB 784748 D02 Section 4.2.

The regulatory information will be displayed in a clear and legible manner without the aid of magnification per KDB 784748 D02 Section 4.3.

The regulatory information is a dedicated exclusion script which is secured and implemented in the factory in an unalterable format. The regulatory information is capable of being retrieved and displayed locally under the full control of the end user. The script is secured to ensure the regulatory information is always displayed correctly.

The following is a sample of the physical label adhered to the exterior of the packaging.

Labels

Model: INX-620

Contains FCC ID: WUW-LWBPLUS

Contains IC: 9613A-LWBPLUS

HVIN: Sterling LWB+

Model: INX-640

Contains FCC ID: WUW-LWBPLUS

Contains IC: 9613A-LWBPLUS

HVIN: Sterling LWB+

Model: INX-660

Contains FCC ID: WUW-LWBPLUS

Contains IC: 9613A-LWBPLUS

HVIN: Sterling LWB+