

Label 1

Label 2


Label 3



# **WARNING** RISK OF ELECTRIC SHOCK



- TO REDUCE THE RISK OF ELECTRIC SHOCK, DO NOT REMOVE FRONT OR BACK COVER.
- READ THE MANUAL AND SAFETY PRECAUTIONS BEFORE PERFORMING INSTALLATION.
- MAKE SURE THE POWER CORD IS PLUGGED INTO THE SYSTEM BEFORE CONNECTING THE POWER CORD TO AN AC POWER OUTLET.
- DISCONNECT POWER FROM THE SYSTEM BY UNPLUGGING THE POWER CORD FROM AC OUTLET.
- THE AC POWER OUTLET SHOULD BE COMPATIBLE WITH THE 3-PIN POWER CORD CONNECTOR.



**Model Name: AID-240SBW**

**Input Rating: 100-240Vac, 2.0-1.0A, 50-60Hz**

For use only with power adapter:

Mfr.: EDAC POWER ELECTRONICS CO LTD

Model: EM11013CWWWW

Input: 100-240Vac, 50-60 Hz, 2.0-1.0A

Output: 19Vdc, 4.73A, Max. 90W

Manufacturer : Avalue Technology Inc.

TEL: +886-2-8226-2345

7F, No.228, Lian Cheng Road, Zhonghe District, New Taipei City, 235, Taiwan



MEDICAL - GENERAL MEDICAL EQUIPMENT AS TO ELECTRICAL SHOCK, FIRE AND MECHANICAL HAZARDS ONLY IN ACCORDANCE WITH ANSI/AAMI ES 60601-1 (2005) + AMD (2012) and CAN/CSA-C22.2 No. 60601-1 (2008) + (2014)

This device complies with Part 15 of the FCC rules. Operation is subject to the following conditions: (1) this device may not cause harmful interference and (2) this device must accept any interference received, including interference that may cause undesired operation



**AID-240SBW Medical Computers**  
FCC ID: XBG-AID240SBWCT IC: 21285-CTNFC  
Contains:



MADE IN TAIWAN

FCC ID: TX2-RTL8821AE IC: 6317A-RTL8821AE

Manufacture Date: 2016-05

Rev.1.30 30ACH

P/N: AID-240SBW01-D02R



S/N: 0000AF00278