



Label Location

Label:

Please adjust it to the same level with the heart. Soft pillow	12.5cm~21.5cm(5.0in~8.4in) 1cm~1.5cm (0.4in~0.6in) Please keep sitting still for 5 mins before starting to measure. This device complies with Part 15 of the FCC Rules. Operation is subject to the following two conditions: (1) This device may not cause harmful interference, and (2) This device must accept any interference received, including interference that may cause undesired operation.
yuwell FCC ID: 2A2JJYE8100B Electronic Blood Pressure Monitor # YE8100B REF YE8100BWH-NENB Power Supply: battery: 3V (2x1.5V AAA) Service life: 5 years (6 times for each day) Made in China JIANGSU YUYUE MEDICAL EQUIPMENT & SUPPLY CO., LTD. No.1 Baisheng Road Development Zone, Danyang, Jiangsu 212300 CHINA	

Note: This label will be permanently engraved on this device.