






Differences Declaration

To whom concern,

We company Xuzhou Yongkang Electronic Science Technology Co.,Ltd. Here by declares:

Model	PI	Bluetooth	auto power off	PCB	Size
YK-80A	Optional	Optional	YES		58*36*32mm
YK-80B	Optional	Optional	—		
YK-80C	Optional	Optional	YES		
YK-81A	Optional	Optional	YES		60*33.5*29mm
YK-81B	Optional	Optional	—		
YK-81C	Optional	Optional	YES		

Model	PI	Bluetooth	auto power off	PCB	Size
YK-82A	Optional	Optional	YES		57.7*35.9*30mm
YK-82B	Optional	Optional	—		
YK-82C	Optional	Optional	YES		
YK-83A	Optional	Optional	YES		58*35.4*31.5mm
YK-83B	Optional	Optional	—		
YK-83C	Optional	Optional	YES		
YK-84A	Optional	Optional	YES		57.9*35.8*30mm
YK-84B	Optional	Optional	—		
YK-84C	Optional	Optional	YES		

1. Description of differences:

The model of Fingertip Pulse Oximeter is: YK-80A, YK-80B, YK-80C, YK-81A, YK-81B, YK-81C, YK-82A, YK-82B, YK-82C, YK-83A, YK-83B, YK-83C, YK-84A, YK-84B, YK-84C, the colors and appearances are different, the others are the same.

2. All of the models have


- (1) Dual color OLED display function
- (2) Blood oxygen measurement and display
- (3) Measure and display pulse rate function
- (4) Low power prompt function

3. POWER

YK-80 series, YK-81 series, YK-82 series, YK-83 series and YK-84 series use two AAA 1.5V alkaline batteries.

4. YK-80 series, YK-81 series, YK-82 series, YK-83 series and YK-84 series use the same PCB board .

Sincerely yours,



Name (in print): Qiu Zhaohao

Company name: Xuzhou Yongkang Electronic Science Technology Co.,Ltd.

Company address: LIANDO U Valley, No.6 Leye Road ,Xuzhou ETDZ, CHINA

Tel: 15062149600

E-Mail: wanxiao@yonker.cn