


Therapy Holdings, Inc.

kGoal Classic / kGoal IntimFlex L3-010

ID Label Information and Location (FCC)

The equipment “can only operate in conjunction with another associated product that has an electronic display” (it is controlled via an app on a smartphone and has no user controls on the device itself). Per FCC and IC guidance, e-labeling is valid for the product.

Please see below screenshot which shows where the app will list FCC ID and model number. The settings page, where the information is listed, is reachable within 1 click of the app’s Home / main page.



For software security, the FCC ID text is stored in the app (hard coded and cannot be modified or overwritten without publishing a new version of the app through the Apple App Store or Google Play Market).

The below figure shows the physical temporary ID label for use on the product packaging per KDB 784748 D02 clause 4.5.

