

1. PLACE ITEM IN POUCH.  
2. PEEL OFF PROTECTIVE STRIP FROM ADHESIVE BAND.  
3. FOLD OVER PAPER FLAP ALONG PERFORATION AND  
PRESS DOWN FIRMLY FROM THE CENTRE OUTWARDS.

Sterilization Date:

1 / 1

Reorder # P23-0149

Stellar Telemetry

Model # 43001-IMP-01  
FCC ID: 2AC4C-AU430001LG  
IC: 12302A-AU430001LG

FCC



Made in USA

Don't use if package is damaged. ISO11607-1 EN868-5 (CE) ②



31/4"x61/2" (83mmx165mm)

Processed yellow

Processed dark grey