

January 19, 2016

Equipment Authorization Division
Federal Communications Commission
7435 Oakland Mills Road
Columbia, MD 21046

FCC ID: 2AD7T11115110401

Product Name: Sense Tag Tag-J120, Tag-J130, Tag-J131

Model difference letter

Models Tag-J120 and Tag-J130 are electrically identical. The only difference is the enclosure. Tag-J120 fits to a standard hospital wristband. Tag-J130 is intended to be used on a lanyard.

Tag-J131 differs electrically only slightly from the two other models, it has cut band tamper detection to detect if a wristband is forcefully removed or cut. The enclosure is designed to fit on a Sonitor wristband.

To locate the electrical difference, please refer to schematic file named "113105-01.Tag-J PBA W LF v3.20.pdf", page 8.

The main PCB in all models is the same.

Sincerely,



January 19, 2016

Hanna Luukkonen
Senior QA and Test Engineer

Date



On behalf of Sonitor Technologies AS