



INTERNAL PHOTOS INF-C396, INF-C398, INF-C399




Figure 1: Photo of the PBA of INF-C39X, front. All the four ultrasound transducer subassemblies seen beside the PBA (one flipped) get mounted to the board for all models. A turret like the one seen to the upper right is mounted to each transducer for INF-C398 and INF-C399 whereas a turret like the one seen to the lower right is mounted to each transducer for INF-C396.




Figure 2: Photo of the PBA of INF-C39X, back.




Figure 3: Photo of the PBA of INF-C39X, with ultrasound transducers mounted, mounted inside the housing.