

JDI 東莞藍創捷特佳電子有限公司

Packing or Production
Instruction

Files No.: RPI-486-XXXX

Page of: 1 0F 2


Date : 2016.5.30

Project No. : 486-15301AN 486-15302AN Product name: GLOW CAP SINGLE ASSEMBLY (BLUE/ORANGE), 100-GCB/O -PDL,RoHS Customer name: Nanthealth

To : R&D/PIE/QA/PROD/AE/QC/IQC/DA

Remarks:

Glow Cap Mac address and SN label Printing content
The labels(2pcs) will be stuck on the main body(1pc)
and be stuck on the specified position of giftbox(1pc)



5. YY means production week
(for example:forty-eighth week
manufacturing is "48")

6. XX means production "year"
(for example :2016 is"16")

1) Nine digit serial number

2) MAC Address(Provided by NH)

3) FCC ID: 2AH33373100 (Provided by NH)

4) Barcode Symbology: QR Code



The SN and MAC
address label will be
stuck here.

Remarks:

- 1.Label P/N:289-51451XN (size:38x12mm)
- 2.The printing Label contents need to be approved by RND.

Efficient

date:: IMMEDIATELY

ZONE	VER	DESCREPTION	DATE	ISSUE BY
	A	New files release	2016.5.30	JM

Issue: JM

Approved by : _____

GlowCap®



NantHealth®
GlowCap®

- + Flashing lights and audible alerts from the GlowCap® provide gentle reminders to take your medication.
- + The GlowCap can be configured to send notifications and reports to family or caregivers to help support your medication adherence.
- + Built-in refill button connects you to your local pharmacy for refill requests.
- + NantOS™ compatible

GlowCap®
100-GCO-PDL
www.nanthhealth.com

Place SN label here

Made in China

NantHealth®

NantHealth®

The SN and MAC
label will be stuck
here