



Antenna Location



Test Positions

Front-of-Face (Separation distance of 10mm)	Back-of-Wristbands (Separation distance of 0mm)
 A photograph showing the front-of-face test position. A wristband is placed on a white foam block, and a red double-headed arrow indicates a separation distance of 10mm between the device and the surface below it.	 A photograph showing the back-of-wristbands test position. A wristband is placed on a white foam block, and a red double-headed arrow indicates a separation distance of 0mm between the device and the surface below it.

Liquid depth

Body Position with 2450MHz liquid depth (15.2cm)	N/A
	N/A