





Fastener location



NOTES:

1. Approximate label size 2.55" x 1.5"
2. Label covers mechanical fasteners



CyMedica Orthopedics® e-vive™ Neuromuscular Electrical Stimulation (NMES) System

Stimulator with protective cover is IP22

Not made with natural rubber latex

Transport and Storage:

-13 °F
(-25 °C)

158 °F
(70 °C)

Temperature:
-13°F to 158°F
(-25°C to 70°C)

10%

90%

Humidity:
10-90%

50kPa

106kPa

Pressure:
50-106kPa

Emergo Europe
Prinsessegracht 20
2514 AP The Hague
The Netherlands

CY-0600-300 Rev. A

Caution: Federal law restricts this device
to sale by or on the order of a physician.

Attention: La loi fédérale restreint vente de
cet appareil par ou sur l'ordre d'un médecin

FCC ID: 2AU28 - CY1000330

Model: e-vive II

This device complies with part 15 of the FCC
Rules. Operation is subject to the following two
conditions: (1) This device may not cause
harmful interference, and (2) this device must
accept any interference received, including
interference that may cause undesired
operation.

Lead wires comply with the
Performance Standard for
electrode lead wires (21
CFR part 898)

ETL CLASSIFIED

Intertek
5000542

Medical Electrical Equipment and Devices

CyMedica Orthopedics
19120 North Pima Road
Suite 135
Scottsdale, AZ 85255 USA
1-844-CYM-2014

1. Approximate label size 2.375" x 5.125"