

M400 Internal Pictures

Cardiac Insight
2375 130th Ave. NE Suite 101
Bellevue, WA USA
FCC ID: 2BEOW-M400

January 30, 2024

Model M400 Sensor • Sensor Module removed from Patient Adhesive Enclosure	M400 Sensor Module  A top-down view of the M400 Sensor Module. The device is oval-shaped with a clear plastic casing. A white label on the left side contains a barcode, the text "CARDEA DUET™", and "Cardiac Insight, Inc. Bellevue, WA 98005 U.S.A.". A ruler is placed below the device, showing it is approximately 3.5 cm long.
Model M400 Sensor • Sensor Module removed from Patient Adhesive Enclosure	Side View  A side view of the M400 Sensor Module, showing its thickness. The device is shown next to a ruler, indicating it is approximately 0.5 cm thick.

Model M400 Sensor • Bottom View, top cover removed. • Gold Contacts are ECG Sensor inputs. • Edge Connector <i>Note: Download connector for ECG data download. When connected to PC the Sensor MCU is powered down, the BT radio is not active.</i>	Bottom View 
Model M400 Sensor • Top View, top cover removed. <i>Note: Download connector for ECG data download. When connected to PC the Sensor MCU is powered down, the BT radio is not active.</i>	Top View 
Model M400 Sensor PCBA • CC2642 BT transceiver • Antenna Match Components • PCB Antenna (Serial Number label removed)	Top Side – Zoom In 

