

Proposed Labeling

This proposed labeling section includes representative copies of the draft product label, draft Instructions for Use, and draft box labeling.

Draft Product Label


Draft product label shown below, to be placed on the back of the C300 Monitor. Red box indicates placement of on-device label (further specified below). Smallest label font size: 4pt.



Product labeling. Graphics color: Pantone Cool Gray 11c, matte finish. Red square indicates placement of on-device label.

Draft On-Device Labeling

On-device label to be positioned on the back of the C300 Monitor, just underneath the USB-C connection port as outlined in the draft product label (above; red rectangle). It will be printed in black on white on a laminated material. The format of the label is as outlined below:




The label will contain:


- GS1 datamatrix, contain all relevant fields (SN#, GTIN, DoM, LOT#)
- SN: 12-character serial number
- GTIN: 14-character GS1 GTIN for C300
- Date of manufacture (MM/DD/YYYY)
- LOT: 4-character LOT number

Power Adapter Labeling

Power adapter uses the original manufacturer labeling with additional Circadia accessory GTIN and LOT #.



Original manufacturer label on adapter



Additional Circadia label on adapter


Draft Box Labeling

Master cartons

The default packaging will be master cartons, in three types of quantities: 10x, 20x and 40x C300 Monitors.

Each carton will have the following external labeling (dimensions in the example below are for the 40 pcs carton):


- 2 long sides will contain Circadia logo + text
- 1st short side will contain Circadia logo only
- 2nd short side will contain regulatory markings, manufacturer information and contents



Draft box labeling, for the 40-pcs carton containing 40 C300 Monitors. Box dimensions and label layout are shown (top), and detailed design of box panels (bottom).

Master Cartons Label - Detailed View

Detailed view of box labeling - 2nd short side of master carton:



The master carton label contains the following (from top to bottom):


- Regulatory symbols
- Manufacturer details
- GTIN barcode for given size of master carton (different between sizes)
- LOT#
- Product reference
- Quantity of what's inside the box

Ziploc Bags


Individual C300 Monitors will be packed in bubble-padded Ziploc bags. Ziploc bags will have individual draft labeling as below:

Front view

Back view



Detailed Ziploc bag labeling:



Contains:
FCC ID: 2BPIS-C300



Manufactured by:
STRIDE Electronics Sdn. Bhd.
No. 50, Lot 1229, Lorong
Perusahaan 2a, Kawasan
Perusahaan Kulim, 09000
Kulim, Kedah,
Malaysia

Manufactured for:
Circadia Health, Inc.
507 S Douglas St.
El Segundo, CA 90245 USA
info@circadia.health
+1 800 985 5596
<https://circadia.health>

Rx Only



REF C300

QTY 1 x Circadia Contactless
Cardiorespiratory Monitor

Designed in California. Assembled in Malaysia.

The position indicated by a red box above in the above label draft will contain a label as below, with fields for:

- GTIN (01)
- Date of manufacture (11)
- Serial number (21)
- LOT number (10)

