



HCT Co., Ltd.
74, Seoicheon-ro 578beon-gil, Majang-myeon, Icheon-si, Gyeonggi-do, 17383
KOREA Tel. +82 31 634 6300 Fax. +82 31 645 6401

FCC ID: A3LSMN981B1

- WPC Test SETUP PHOTOGRAPHS

REVISION HISTORY


The revision history for this document is shown in table.

Revision No.	Date of Issue	Description
0	Oct 29, 2021	Initial Release
1	Nov.11,2021	Revised the section description.

1. Test setup Photo of Direct contact from the edge of the Device to the center of the Probe :5.4 cm
Rear / Front side



Left / Right side



Top/ Bottom side




2. Test setup Photo of Distance from the rear side of the Device to the center of the Probe from 6.4 cm to 15.4 cm

6.4cm,7.4cm

Distance from the rear side of the Device to the center of the Probe from 6.4 cm



Distance from the rear side of the Device to the center of the Probe from 7.4 cm




8.4cm. 9.4cm

Distance from the rear side of the Device to the center of the Probe from 8.4 cm



Distance from the rear side of the Device to the center of the Probe from 9.4 cm




10.4cm. 11.4cm

Distance from the rear side of the Device to the center of the Probe from 10.4 cm




Distance from the rear side of the Device to the center of the Probe from 11.4 cm



12.4cm. 13.4cm

Distance from the rear side of the Device to the center of the Probe from 12.4 cm




Distance from the rear side of the Device to the center of the Probe from 13.4 cm



14.4cm. 15.4cm

Distance from the rear side of the Device to the center of the Probe from 14.4 cm



Distance from the rear side of the Device to the center of the Probe from 15.4 cm

