

Appendix to Report

Dosimetric Assessment of the Portable Device Panasonic KX-TGA20 (FCC ID: ACJ96NKX-TGA20) (IC 216A-KXTGA20)

According to the FCC Requirements

December 18, 2012

IMST GmbH
Carl-Friedrich-Gauß-Str. 2
D-47475 Kamp-Lintfort

Customer
Panasonic System networks Co., Ltd.
1-62, 4-come, Minoshima, Hakata-ku
5215 MC's-Hertogenbosch
Fukuoka 812-5831
Japan

The test results only relate to the items tested. This report shall not be reproduced except in full without the written approval of the testing laboratory.

1 Subject of Investigation



Fig. 1: Picture of the device under test.

8.9 Pictures of the Device under Test

Figure 12 shows the device under test.



Fig. 12: Picture of the device under test.

8.10 Test Positions for the Device under Test

Figure 13 - 14 show the test positions for the SAR measurements.



Fig. 13: Body worn configuration for the Panasonic KX-TGA20, front of the device towards the phantom.



Fig. 14: Body worn configuration for the Panasonic KX-TGA20, back of the device towards the phantom.