

IDENTIFICATION LABEL:

PROPOSED LABEL

The shall be permanently affixed in the rear side of the equipment

BP170B(Blood Pressure Monitor)      Input Voltage: DC 6V$\pm$, 2A or 1.5V AA Battery 4EA Measuring Range: 0~300mmHg InBody  InBody Co., Ltd. 15, Heugam-gil, Ipjang-myeon, Seobuk-gu, Cheonan-si, Chungcheongnam-do 31025 KOREA TEL: +82-41-581-3003 FAX: +82-41-581-3103 Website: inbody.com E-mail: info@inbody.com [EC REP] InBody Europe B.V. Gyroscopweg 122, 1042 AZ, Amsterdam, The Netherlands TEL: +31-20-238-6080 FAX: +31-6-5734-1858 Website: nl.inbody.com E-mail: info.eu@inbody.com L5-F-109-E-200222 MADE IN KOREA FCC ID: F6O-INBODY-BP170B IC: 22967-BP170B HVIN: BP170B	BP160B(Blood Pressure Monitor)      Input Voltage: DC 6V$\pm$, 2A or 1.5V AA Battery 4EA Measuring Range: 0~300mmHg InBody  InBody Co., Ltd. 15, Heugam-gil, Ipjang-myeon, Seobuk-gu, Cheonan-si, Chungcheongnam-do 31025 KOREA TEL: +82-41-581-3003 FAX: +82-41-581-3103 Website: inbody.com E-mail: info@inbody.com [EC REP] InBody Europe B.V. Gyroscopweg 122, 1042 AZ, Amsterdam, The Netherlands TEL: +31-20-238-6080 FAX: +31-6-5734-1858 Website: nl.inbody.com E-mail: info.eu@inbody.com L5-F-109-E-200222 MADE IN KOREA FCC ID: F6O-INBODY-BP170B IC: 22967-BP170B HVIN: BP160B
--	--

