

1. THIS SPECIFICATION IS CONTROLLED BY CRHF PACKAGE ENGINEERING. THIS SPECIFICATION IS INTENDED TO BE THE FINAL LABEL DESIGN OUTPUT GENERATED FROM A COMBINATION OF LABEL PARTS WITHIN THE ASSOCIATED PDM BOM.
2. THESE LABEL IMAGES ARE THE COMBINATION OF THE TEMPLATE FILE, MODEL DATA RECORD, SAMPLE TRIGGER DATA AND REPRESENTATIVE VIEWS OF THE LABEL DIE LINES.
3. PRINTING REQUIREMENTS - THE LABELS ARE TO BE PRINTED USING A 300 DPI ZEBRA THERMAL TRANSFER PRINTER.
4. THE LOT OR SERIAL NUMBER, USE BY DATE, DATE OF MANUFACTURE, PIN NUMBER, ADDRESS CODE, GTIN, PRODUCT NUMBER (REF) AND BAR CODES ARE SHOWN FOR EXAMPLE ONLY. THESE FIELDS ARE VARIABLE AND ARE TO BE SUPPLIED BY MANUFACTURING AT THE TIME OF PRINTING.
5. SPACING AND TEXT WRAPPING SHOWN IN THE IMAGES MAY DIFFER SLIGHTLY FROM ACTUAL PRINTED LABEL, AS LONG AS CONTENT OF PRINTED VERSION IS ACCURATE. THE IMAGES INCLUDE REPRESENTATION OF THE PRE-PRINTED LABEL STOCK. BRANDING, PART NUMBER BAR CODES (2D OR OTHER) AND OTHER PRE-PRINTED AREAS OR ADDITIONAL STICKERS MAY CHANGE AND ARE ACCEPTABLE.
6. THE LABEL PRINTED IN MANUFACTURING MUST BE AT AN EQUAL OR LATER MODEL DATA AND TEMPLATE REVISION AS SHOWN ON THE IMAGE.
7. (FOR LABEL SET ONLY) THE TRIGGER DATA SHOWN ON THE LABEL SET IMAGE MAY NOT MATCH THE PRINTED LABEL.

SPEC TAB	MFG ADDRESS
A86277-001	Tolochenaz, Switzerland
A86277-002	Juncos, Puerto Rico, USA

DeviceImageSpec_02

Medtronic MEDTRONIC, INC. MINNEAPOLIS, MINNESOTA USA MEDTRONIC CONFIDENTIAL THIS DOCUMENT IS THE PROPERTY OF MEDTRONIC, INC. AND MUST BE ACCOUNTED FOR. INFORMATION HEREON IS CONFIDENTIAL. DO NOT REPRODUCE IT, REVEAL IT TO UNAUTHORIZED PERSONS OR SEND IT OUTSIDE MEDTRONIC WITHOUT PROPER AUTHORIZATION.	TITLE LABEL IMAGE SPEC-SERENA QUAD CRT-P MRI SURESCAN W4TR02 US				
	LABEL IMAGE SPEC TEMPLATE: A83320 REV C				
	DATE 07-Feb-2017	SCALE N/A	ENGINEERING JUDY SALZER		
	MODEL W4TR02	Size B	PART NUMBER A86277	REV B	SHEET 1 of 3

Serena™ Quad CRT-P MRI SureScan™

Bluetooth™

OAE-DDDR

ACCELEROMETER

MR Conditional

CRT-P

Model **W4TR02**

SN WWW000000W

2000-00-00
YYYY-MM-DD

2000-00-00
YYYY-MM-DD

REF W4TR02

Package contents: 1x

Implantable device

Torque wrench

Product documentation

DDD

↑ 130 bpm
↓ 50 bpm

Upper tracking rate and lower rate

0.3 mV

Sensitivity: RA

0.9 mV

Sensitivity: RV

→ 130 ms
→ 100 ms

A-V interval (paced and sensed)

0.4 ms

Amplitude and pulse width: RA

3.5 V
→ 0.4 ms

Amplitude and pulse width: RV

3.5 V
→ 0.4 ms

Amplitude and pulse width: LV

LV1 - LV2

Pacing polarity: LV

LV->RV
0 ms

V pacing and V-V pace delay

AUTO

Pacing polarity: RA, RV

Adaptive

AUTO

Sensing polarity: RA, RV

OFF

Therapies: AT, AF

AUTO

PVAPR

Rx Only

FCFID: LFSBLEIMPLANT2

Medtronic, Inc., 710 Medtronic Parkway, Minneapolis, MN 55432 USA

Manufactured in: Tolochenaz, Switzerland

www.medtronic.com

STERILE EO

18 °C
0 °F

55 °C
131 °F

i

Serena™ Quad CRT-P MRI SureScan™

Bluetooth™

RA, IS-1
RV, IS-1
OAE-DDDR
LV, IS4 - LLLL
ACCELEROMETER

CRT-P

W4TR02
SN WWW000000W
2000-00-00
YYYY-MM-DD
2000-00-00
YYYY-MM-DD
REF W4TR02

Burst+/- RA
Ramp: RA

Package contents: 1x

Implantable device
Torque wrench

DDD

↑ 130 bpm
↓ 50 bpm
Upper tracking rate and lower rate

0.3 mV
Sensitivity: RA

0.9 mV
Sensitivity: RV

→ 130 ms
→ 100 ms
A-V interval (paced and sensed)

↑ 3.5 V
↔ 0.4 ms
Amplitude and pulse width: RA

↑ 3.5 V
↔ 0.4 ms
Amplitude and pulse width: RV

↑ 3.5 V
↔ 0.4 ms
Amplitude and pulse width: LV

LV1 - LV2
Pacing polarity: LV

LV->RV
0 ms
V pacing and V-V pace delay

AUTO
Pacing polarity: RA, RV

Adaptive

AUTO
Sensing polarity: RA, RV

OFF
Therapies: AT, AF

AUTO
PVARP

Rx Only

FCC ID: LFSBLEIMPLANT2

STERILE EO

55 °C
131 °F

-18 °C
0 °F


Medtronic, Inc., 710 Medtronic Parkway, Minneapolis, MN 55432 USA
Manufactured in: Tolochenaz, Switzerland
www.medtronic.com

M968371A001 B - M970302A001 A

M968371A001 B - M970302A001 A		
		2000-00-00 YYYY-MM-DD
Medtronic W4TR02	 SureScan	 Serena™ Quad CRT-P MRI SureScan™ SN WWW000000W
Medtronic W4TR02	 SureScan	 Serena™ Quad CRT-P MRI SureScan™ SN WWW000000W
Medtronic W4TR02	 SureScan	 Serena™ Quad CRT-P MRI SureScan™ SN WWW000000W
Medtronic W4TR02	 SureScan	 Serena™ Quad CRT-P MRI SureScan™ SN WWW000000W
Medtronic W4TR02	 SureScan	 Serena™ Quad CRT-P MRI SureScan™ SN WWW000000W
Medtronic W4TR02	 SureScan	 Serena™ Quad CRT-P MRI SureScan™ SN WWW000000W
Medtronic W4TR02	 SureScan	 Serena™ Quad CRT-P MRI SureScan™ SN WWW000000W
Medtronic W4TR02	 SureScan	 Serena™ Quad CRT-P MRI SureScan™ SN WWW000000W
Medtronic W4TR02	 SureScan	 Serena™ Quad CRT-P MRI SureScan™ SN WWW000000W
Medtronic W4TR02	 SureScan	 Serena™ Quad CRT-P MRI SureScan™ SN WWW000000W
Medtronic W4TR02	 SureScan	 Serena™ Quad CRT-P MRI SureScan™ SN WWW000000W
Medtronic W4TR02	 SureScan	 Serena™ Quad CRT-P MRI SureScan™ SN WWW000000W

Medtronic Serena™ Quad CRT-P MRI SureScan™ W4TR02 Implantation date:		 Bluetooth®
	2000-00-00 YYYY-MM-DD	 SureScan
 MR Conditional	REF W4TR02 	
SN WWW000000W 		

Medtronic Serena™ Quad CRT-P MRI SureScan™ W4TR02 Implantation date:		 Bluetooth®
	2000-00-00 YYYY-MM-DD	 SureScan
 MR	REF W4TR02 	
SN WWW000000W 		

Medtronic Serena™ Quad CRT-P MRI SureScan™ W4TR02 Implantation date:		 Bluetooth®
	2000-00-00 YYYY-MM-DD	 SureScan
 MR	REF W4TR02 	
SN WWW000000W 		








Serena™ Quad CRT-P MRI SureScan™ W4TR02

- AdaptivCRT™ Algorithm
- VectorExpress™ LV Automated Test
- PhysioCurve™ Design
- OptiVol™ 2.0 Fluid Status Monitoring

REF W4TR02
SN WWW000000W
2000-00-00
YYYY-MM-DD

M968371A001 B - M970302A001 A

SPEC TAB
A86277-001

NOTE: IMAGES NOT TO SCALE
 LABEL IMAGE SPEC TEMPLATE: A83320 REV C

PART NUMBER	REV	SHEET
A86277	B	2 of 3



Medtronic


Bluetooth

Serena™ Quad CRT-P MRI
SureScan™



CRT-P



RA, IS-1
RV, IS-1
LV, IS4 - LLLL
OAE-DDDR
ACCELEROMETER
SureScan
MR Conditional

(Model)	W4TR02
SN	WWW000000W
	2000-00-00 YYYY-MM-DD
	2000-00-00 YYYY-MM-DD
REF	W4TR02

Burst+: RA
Ramp: RA

Package contents: 1x
Implantable device
Torque wrench
Product documentation

DDD



↑ 130 bpm
↓ 50 bpm
Upper tracking rate and lower rate



0.3 mV
Sensitivity: RA



0.9 mV
Sensitivity: RV



→ 130 ms
→ 100 ms
A-V interval (paced and sensed)



↑ 3.5 V
↔ 0.4 ms
Amplitude and pulse width: RA



↑ 3.5 V
↔ 0.4 ms
Amplitude and pulse width: RV



↑ 3.5 V
↔ 0.4 ms
Amplitude and pulse width: LV



LV1 - LV2
Pacing polarity: LV



LV->RV
0 ms
V pacing and V-V pace delay



AUTO
Pacing polarity: RA, RV

Adaptive



AUTO
Sensing polarity: RA, RV
Therapies: AT, AF
PVARP

Rx Only

FCFID: LFSBLEIMPLANT2

STERILE EO



55 °C
131 °F
-18 °C
0 °F



Medtronic, Inc., 710 Medtronic Parkway, Minneapolis, MN 55432 USA
Manufactured in: Juncos, Puerto Rico, USA
Product of USA
www.medtronic.com

M968371A001 B - M970302A001 A






Medtronic


Bluetooth™

Serena™ Quad CRT-P MRI SureScan™



RA, IS-1
RV, IS-1
LV, IS4 - LLLL

OAE-DDDR



ACCELEROMETER



SureScan
MR Conditional



Burst+/- RA



Ramp: RA



CRT-P



Package contents: 1x



Implantable device



Torque wrench

DDD



↑ 130 bpm
↓ 50 bpm

Upper tracking rate and lower rate



0.3 mV

Sensitivity: RA



0.9 mV

Sensitivity: RV



→ 130 ms
→ 100 ms

A-V interval (paced and sensed)



↑ 3.5 V
↔ 0.4 ms

Amplitude and pulse width: RA



↑ 3.5 V
↔ 0.4 ms

Amplitude and pulse width: RV



↑ 3.5 V
↔ 0.4 ms

Amplitude and pulse width: LV



LV1 - LV2

Pacing polarity: LV



LV->RV
0 ms

V pacing and V-V pace delay



AUTO

Pacing polarity: RA, RV



Adaptive



AUTO

Sensing polarity: RA, RV



OFF

Therapies: AT, AF



AUTO

PVARP

Rx Only

FCC ID: LFSBLEIMPLANT2

STERILE

EO




-18 °C
0 °F

55 °C
131 °F



Medtronic, Inc., 710 Medtronic Parkway, Minneapolis, MN 55432 USA

Manufactured in: Juncos, Puerto Rico, USA

www.medtronic.com




M968371A001 B - M970302A001 A

M968371A001 B - M970302A001 A		
		2000-00-00 YYYY-MM-DD
Medtronic W4TR02	 SureScan	 Serena™ Quad CRT-P MRI SureScan™ SN WWW000000W
Medtronic W4TR02	 SureScan	 Serena™ Quad CRT-P MRI SureScan™ SN WWW000000W
Medtronic W4TR02	 SureScan	 Serena™ Quad CRT-P MRI SureScan™ SN WWW000000W
Medtronic W4TR02	 SureScan	 Serena™ Quad CRT-P MRI SureScan™ SN WWW000000W
Medtronic W4TR02	 SureScan	 Serena™ Quad CRT-P MRI SureScan™ SN WWW000000W
Medtronic W4TR02	 SureScan	 Serena™ Quad CRT-P MRI SureScan™ SN WWW000000W
Medtronic W4TR02	 SureScan	 Serena™ Quad CRT-P MRI SureScan™ SN WWW000000W
Medtronic W4TR02	 SureScan	 Serena™ Quad CRT-P MRI SureScan™ SN WWW000000W
Medtronic W4TR02	 SureScan	 Serena™ Quad CRT-P MRI SureScan™ SN WWW000000W
Medtronic W4TR02	 SureScan	 Serena™ Quad CRT-P MRI SureScan™ SN WWW000000W
Medtronic W4TR02	 SureScan	 Serena™ Quad CRT-P MRI SureScan™ SN WWW000000W
Medtronic W4TR02	 SureScan	 Serena™ Quad CRT-P MRI SureScan™ SN WWW000000W

Medtronic Serena™ Quad CRT-P MRI SureScan™ W4TR02 Implantation date:		 Bluetooth®
	2000-00-00 YYYY-MM-DD	 SureScan
 MR Conditional	REF W4TR02 	
SN WWW000000W 		

Medtronic Serena™ Quad CRT-P MRI SureScan™ W4TR02 Implantation date:		 Bluetooth®
	2000-00-00 YYYY-MM-DD	 SureScan
 MR	REF W4TR02 	
SN WWW000000W 		

Medtronic Serena™ Quad CRT-P MRI SureScan™ W4TR02 Implantation date:		 Bluetooth®
	2000-00-00 YYYY-MM-DD	 SureScan
 MR	REF W4TR02 	
SN WWW000000W 		








Serena™ Quad CRT-P MRI SureScan™ W4TR02

- AdaptivCRT™ Algorithm
- VectorExpress™ LV Automated Test
- PhysioCurve™ Design
- OptiVol™ 2.0 Fluid Status Monitoring

REF W4TR02
SN WWW000000W
2000-00-00
YYYY-MM-DD

M968371A001 B - M970302A001 A

SPEC TAB
A86277-002

NOTE: IMAGES NOT TO SCALE
 LABEL IMAGE SPEC TEMPLATE: A83320 REV C

PART NUMBER	REV	SHEET
A86277	B	3 of 3