

1. THIS SPECIFICATION IS CONTROLLED BY CRHF PACKAGE ENGINEERING. THIS SPECIFICATION IS INTENDED TO BE THE FINAL LABEL DESIGN OUTPUT GENERATED FROM A COMBINATION OF LABEL PARTS WITHIN THE ASSOCIATED PDM BOM.
2. THESE LABEL IMAGES ARE THE COMBINATION OF THE TEMPLATE FILE, MODEL DATA RECORD, SAMPLE TRIGGER DATA AND REPRESENTATIVE VIEWS OF THE LABEL DIE LINES.
3. PRINTING REQUIREMENTS - THE LABELS ARE TO BE PRINTED USING A 300 DPI ZEBRA THERMAL TRANSFER PRINTER.
4. THE LOT OR SERIAL NUMBER, USE BY DATE, DATE OF MANUFACTURE, PIN NUMBER, ADDRESS CODE, GTIN, PRODUCT NUMBER (REF) AND BAR CODES ARE SHOWN FOR EXAMPLE ONLY. THESE FIELDS ARE VARIABLE AND ARE TO BE SUPPLIED BY MANUFACTURING AT THE TIME OF PRINTING.
5. SPACING AND TEXT WRAPPING SHOWN IN THE IMAGES MAY DIFFER SLIGHTLY FROM ACTUAL PRINTED LABEL, AS LONG AS CONTENT OF PRINTED VERSION IS ACCURATE. THE IMAGES INCLUDE REPRESENTATION OF THE PRE-PRINTED LABEL STOCK. BRANDING, PART NUMBER BAR CODES (2D OR OTHER) AND OTHER PRE-PRINTED AREAS OR ADDITIONAL STICKERS MAY CHANGE AND ARE ACCEPTABLE.
6. THE LABEL PRINTED IN MANUFACTURING MUST BE AT AN EQUAL OR LATER MODEL DATA AND TEMPLATE REVISION AS SHOWN ON THE IMAGE.
7. (FOR LABEL SET ONLY) THE TRIGGER DATA SHOWN ON THE LABEL SET IMAGE MAY NOT MATCH THE PRINTED LABEL.

SPEC TAB	MFG ADDRESS
A86279-001	Tolochenaz, Switzerland
A86279-002	Juncos, Puerto Rico, USA


DeviceImageSpec_02

Medtronic MEDTRONIC, INC. MINNEAPOLIS, MINNESOTA USA MEDTRONIC CONFIDENTIAL THIS DOCUMENT IS THE PROPERTY OF MEDTRONIC, INC. AND MUST BE ACCOUNTED FOR. INFORMATION HEREON IS CONFIDENTIAL. DO NOT REPRODUCE IT, REVEAL IT TO UNAUTHORIZED PERSONS OR SEND IT OUTSIDE MEDTRONIC WITHOUT PROPER AUTHORIZATION.	TITLE LABEL IMAGE SPEC-SOLARA QUAD CRT-P MRI SURESCAN W4TR03 US				
	LABEL IMAGE SPEC TEMPLATE: A83320 REV C				
	DATE 07-Feb-2017	SCALE N/A	ENGINEERING JUDY SALZER		
	MODEL W4TR03	Size B	PART NUMBER A86279	REV B	SHEET 1 of 3



Solara™ Quad CRT-P MRI SureScan™










CRT-P



MR Conditional






Model **W4TR03**

SN WWW000000W

2000-00-00
YYYY-MM-DD

2000-00-00
YYYY-MM-DD

REF W4TR03

2000-00-00
YYYY-MM-DD

DDD



↑ 130 bpm
↓ 50 bpm

Upper tracking rate and lower rate



0.3 mV

Sensitivity: RA



0.9 mV

Sensitivity: RV



→ 130 ms
→ 100 ms

A-V interval (paced and sensed)



↑ 3.5 V
↔ 0.4 ms

Amplitude and pulse width: RA



↑ 3.5 V
↔ 0.4 ms

Amplitude and pulse width: RV



↑ 3.5 V
↔ 0.4 ms

Amplitude and pulse width: LV



LV1 - LV2

Pacing polarity: LV



LV->RV
0 ms

V pacing and V-V pace delay



AUTO

Pacing polarity: RA, RV

Adaptive



AUTO

Sensing polarity: RA, RV



OFF

Therapies: AT, AF



AUTO

PVARP

Rx Only

FCFID: LFSBLEIMPLANT2

Medtronic, Inc., 710 Medtronic Parkway, Minneapolis, MN 55432 USA
Manufactured in: Tolochenaz, Switzerland
www.medtronic.com

STERILE EO



55 °C
131 °F



-18 °C
0 °F







Medtronic


Bluetooth™

Solara™ Quad CRT-P


MRI SureScan™



CRT-P



SureScan™



(Mode) **W4TR03**

SN WWW000000W

 **2000-00-00**
YYYY-MM-DD

 2000-00-00
YYYY-MM-DD

REF W4TR03



Burst+/- RA



Ramp: RA



Package contents:

1x



Implantable device



Torque wrench

DDD



↑ 130 bpm
↓ 50 bpm

Upper tracking rate and lower rate



0.3 mV

Sensitivity: RA



0.9 mV

Sensitivity: RV



→ 130 ms
→ 100 ms

A-V interval (paced and sensed)



↑ 3.5 V
↔ 0.4 ms

Amplitude and pulse width: RA



↑ 3.5 V
↔ 0.4 ms

Amplitude and pulse width: RV



↑ 3.5 V
↔ 0.4 ms

Amplitude and pulse width: LV



LV1 - LV2

Pacing polarity: LV



LV->RV
0 ms

V pacing and V-V pace delay



AUTO

Pacing polarity: RA, RV



Adaptive



AUTO

Sensing polarity: RA, RV



OFF

Therapies: AT, AF



AUTO

PVARP

Rx Only

FCC ID: LFS5BLEIMPLANT2

STERILE

EO




-18 °C
0 °F



55 °C
131 °F



Medtronic, Inc., 710 Medtronic Parkway, Minneapolis, MN 55432 USA

Manufactured in: Tolochenaz, Switzerland

www.medtronic.com



M968371A001 B - M970304A001 A

M968371A001 B - M970304A001 A		
		2000-00-00 YYYY-MM-DD
Medtronic W4TR03	 SureScan	 Solara™ Quad CRT-P MRI SureScan™ SN WWW000000W
Medtronic W4TR03	 SureScan	 Solara™ Quad CRT-P MRI SureScan™ SN WWW000000W
Medtronic W4TR03	 SureScan	 Solara™ Quad CRT-P MRI SureScan™ SN WWW000000W
Medtronic W4TR03	 SureScan	 Solara™ Quad CRT-P MRI SureScan™ SN WWW000000W
Medtronic W4TR03	 SureScan	 Solara™ Quad CRT-P MRI SureScan™ SN WWW000000W
Medtronic W4TR03	 SureScan	 Solara™ Quad CRT-P MRI SureScan™ SN WWW000000W
Medtronic W4TR03	 SureScan	 Solara™ Quad CRT-P MRI SureScan™ SN WWW000000W
Medtronic W4TR03	 SureScan	 Solara™ Quad CRT-P MRI SureScan™ SN WWW000000W
Medtronic W4TR03	 SureScan	 Solara™ Quad CRT-P MRI SureScan™ SN WWW000000W
Medtronic W4TR03	 SureScan	 Solara™ Quad CRT-P MRI SureScan™ SN WWW000000W
Medtronic W4TR03	 SureScan	 Solara™ Quad CRT-P MRI SureScan™ SN WWW000000W
Medtronic W4TR03	 SureScan	 Solara™ Quad CRT-P MRI SureScan™ SN WWW000000W

Medtronic

Solara™ Quad CRT-P MRI SureScan™ W4TR03

Implantation date:

2000-00-00

YYYY-MM-DD

REF

W4TR03



SN

WWW000000W



Medtronic

Solara™ Quad CRT-P MRI SureScan™ W4TR03

Implantation date:

2000-00-00

YYYY-MM-DD

SN

WWW000000W



Medtronic

Solara™ Quad CRT-P MRI SureScan™ W4TR03

Implantation date:

2000-00-00

YYYY-MM-DD

SN

WWW000000W



Medtronic

OAE-DDDR

2x IS-1 BI 1x IS4-LLLL

Solara™ Quad CRT-P MRI SureScan™ W4TR03

REF W4TR03

SN WWW000000W

2000-00-00
YYY-MM-DD

(01)00000000000000(17)000000(21)WWW000000W

M968371A001 B - M970304A001 A

End Label 1





OAE-DDDR

2x

IS-1 BI

1x

IS4-LLLL

Solara™ Quad CRT-P MRI SureScan™ W4TR03

- VectorExpress™ LV Automated Test
- PhysioCurve™ Design
- OptiVol™ 2.0 Fluid Status Monitoring

REF W4TR03

SN WWW000000W

2000-00-00

YYYY-MM-DD

MR Conditional

PIN: V0W000000W

End Label 2

SPEC TAB
A86279-001

A86279-001

NOTE: IMAGES NOT TO SCALE
 LABEL IMAGE SPEC TEMPLATE: A83320 REV C


DeviceImageSpec 02

PART NUMBER	REV	SHEET
A86279	B	2 of 3



Solara™ Quad CRT-P MRI SureScan™


















DDD



↑ 130 bpm
↓ 50 bpm

Upper tracking rate and lower rate



0.3 mV

Sensitivity: RA



0.9 mV

Sensitivity: RV



→ 130 ms
→ 100 ms

A-V interval (paced and sensed)



↑ 3.5 V
↔ 0.4 ms

Amplitude and pulse width: RA



↑ 3.5 V
↔ 0.4 ms

Amplitude and pulse width: RV



↑ 3.5 V
↔ 0.4 ms

Amplitude and pulse width: LV



LV1 - LV2

Pacing polarity: LV



LV->RV
0 ms

V pacing and V-V pace delay



AUTO

Pacing polarity: RA, RV



Adaptive



AUTO

Sensing polarity: RA, RV



OFF

Therapies: AT, AF



AUTO

PVARP

Rx Only

FCC ID: LFSBLEIMPLANT2

STERILE

EO





55 °C
131 °F

-18 °C
0 °F



Medtronic, Inc., 710 Medtronic Parkway, Minneapolis, MN 55432 USA
Manufactured in: Juncos, Puerto Rico, USA
www.medtronic.com

Product of USA



M968371A001 B - M970304A001 A



Medtronic


Bluetooth™

Solara™ Quad CRT-P


MRI SureScan™



CRT-P



MRI SureScan



(Mode) **W4TR03**

SN WWW000000W

2000-00-00
YYYY-MM-DD

2000-00-00
YYYY-MM-DD

REF W4TR03



Burst + RA



Ramp: RA



Package contents:

1x



Implantable device



Torque wrench

DDD



↑ 130 bpm
↓ 50 bpm

Upper tracking rate and lower rate



0.3 mV

Sensitivity: RA



0.9 mV

Sensitivity: RV



→ 130 ms
→ 100 ms

A-V interval (paced and sensed)



↑ 3.5 V
↔ 0.4 ms

Amplitude and pulse width: RA



↑ 3.5 V
↔ 0.4 ms

Amplitude and pulse width: RV



↑ 3.5 V
↔ 0.4 ms

Amplitude and pulse width: LV



LV1 - LV2

Pacing polarity: LV



LV->RV
0 ms

V pacing and V-V pace delay



AUTO

Pacing polarity: RA, RV



Adaptive



AUTO

Sensing polarity: RA, RV



OFF

Therapies: AT, AF



AUTO

PVARP

Rx Only

FCC ID: LFSBLEIMPLANT2

STERILE

EO




-18 °C
0 °F



55 °C
131 °F



Medtronic, Inc., 710 Medtronic Parkway, Minneapolis, MN 55432 USA

Manufactured in: Juncos, Puerto Rico, USA

www.medtronic.com

M968371A001 B - M970304A001 A



M968371A001 B - M970304A001 A		
		2000-00-00 YYYY-MM-DD
Medtronic W4TR03	 SureScan	 Solara™ Quad CRT-P MRI SureScan™ SN WWW000000W
Medtronic W4TR03	 SureScan	 Solara™ Quad CRT-P MRI SureScan™ SN WWW000000W
Medtronic W4TR03	 SureScan	 Solara™ Quad CRT-P MRI SureScan™ SN WWW000000W
Medtronic W4TR03	 SureScan	 Solara™ Quad CRT-P MRI SureScan™ SN WWW000000W
Medtronic W4TR03	 SureScan	 Solara™ Quad CRT-P MRI SureScan™ SN WWW000000W
Medtronic W4TR03	 SureScan	 Solara™ Quad CRT-P MRI SureScan™ SN WWW000000W
Medtronic W4TR03	 SureScan	 Solara™ Quad CRT-P MRI SureScan™ SN WWW000000W
Medtronic W4TR03	 SureScan	 Solara™ Quad CRT-P MRI SureScan™ SN WWW000000W
Medtronic W4TR03	 SureScan	 Solara™ Quad CRT-P MRI SureScan™ SN WWW000000W
Medtronic W4TR03	 SureScan	 Solara™ Quad CRT-P MRI SureScan™ SN WWW000000W
Medtronic W4TR03	 SureScan	 Solara™ Quad CRT-P MRI SureScan™ SN WWW000000W
Medtronic W4TR03	 SureScan	 Solara™ Quad CRT-P MRI SureScan™ SN WWW000000W

Medtronic

Solara™ Quad CRT-P MRI SureScan™ W4TR03

Implantation date:

2000-00-00

YYYY-MM-DD

SureScan

MR Conditional

Bluetooth

REF

W4TR03



SN

WWW000000W



Medtronic

Solara™ Quad CRT-P MRI SureScan™ W4TR03

Implantation date:

2000-00-00

YYYY-MM-DD

SureScan

MR

Bluetooth

SN

WWW000000W



Medtronic

Solara™ Quad CRT-P MRI SureScan™ W4TR03

Implantation date:

2000-00-00

YYYY-MM-DD

SureScan

MR

Bluetooth

SN

WWW000000W





Medtronic

Solara™ Quad CRT-P MRI SureScan™ W4TR03

OAE-DDDR

2x IS-1 BI **1x IS4-LLLL**



(01)00000000000000(17)000000(21)WWW000000W

REF W4TR03

SN WWW000000W



2000-00-00
YYYY-MM-DD



M968371A001 B - M970304A001 A



OAE-DDDR

2x

IS-1 BI

1x

IS4-LLLL

Solara™ Quad CRT-P MRI SureScan™ W4TR03

- VectorExpress™ LV Automated Test
- PhysioCurve™ Design
- OptiVol™ 2.0 Fluid Status Monitoring



REF W4TR03

SN WWW000000W



2000-00-00
YYYY-MM-DD

M968371A001 B - M970304A001 A

MR Conditional

PIN: VV0W000000W

SPEC TAB
A86279-002

NOTE: IMAGES NOT TO SCALE
 LABEL IMAGE SPEC TEMPLATE: A83320 REV C

PART NUMBER	REV	SHEET
A86279	B	3 of 3