

1. Test Setup Photo


1. Test Setup Photo



1. Test Setup Photo




1. Test Setup Photo

1. Test Setup Photo



1. Test Setup Photo

Test Setup Photo



Test Setup Photo

[835MHz Head / Body Tissue Liquid Depth : 15cm]




Test Setup Photo**[Conducted Power Measurement set-up]**

Conducted power test

TDC-6610 Conducted Power test

CDMA (363 channel)	
Time	TX Power (dBm)
0	26.00
2	25.98
4	25.96
6	25.94
8	25.91
10	25.88
12	25.85
14	25.82
16	25.80
18	25.78
20	25.76
22	25.74
24	25.72
26	25.71
28	25.68
30	25.66



1. Test Setup Photo (Validation test)