



Internal ref: ICON2/LET/0001/PB

Lund 03/17/2006

This is a description of the major differences between the original product Precise 250MC and the new version Precise 200MC. Both will be on the market targeting different customers.

<u>Part type</u>	<u>Precise 250MC (MS 010 114)</u>	<u>Precise 200MC (MS 010 126)</u>
Sensor	UPEK TCS1CD	UPEK TCS2CF (smaller sensor area)
Front cover	MB 010 047 (plastic)	MB 010 054 (plastic)
Smart card guide	MB 010 048 (plastic)	-
Sensor gasket	MB 010 050 (plastic)	-
Rear cover	MB 010 046 (metal)	MB 010 053 (plastic)

In my opinion, the only difference above that will affect radiated emission is the rear cover. This is a metal part for 250MC and a plastic one for 200MC.

Regards,
Per-Ola Olsson

Phone: +46 46 31 11 00
Fax: +46 46 31 11 01
E-mail: po@precisebiometrics.com
Web: www.precisebiometrics.com

Scheelevägen 19 C
SE-223 70 Lund
Sweden