



F.1 Test Sample



F.2 Test positions

Top Edge	Back Face
 A photograph showing the top edge of a device being tested. The device is mounted on a stand, and a computer monitor displays test results. The setup is in a laboratory environment.	 A photograph showing the back face of a device being tested. The device is mounted on a stand, and a computer monitor displays test results. The setup is in a laboratory environment.

F.3 Antenna host platform location and adjacent edge positions relative to the body



F.4 Phantom liquid level during measurements

ELI phantom Body liquid SAR Dasy SPEAG MBBL600-6000V6 Batch - 160630-01

