



1. Test Setup Photo



Test Setup Photo




Test Setup Photo



[Distance of 25 mm from the EUT]

Test Setup Photo


[835MHz Body Tissue Liquid Depth : 15cm]



! Battery Power test

■ AXW-P800B Conducted Power test

CDMA (363 channel)	
Time	TX Power (dBm)
0	24.50
2	24.49
4	24.47
6	24.45
8	24.42
10	24.41
12	24.40
14	24.39
16	24.38
18	24.37
20	24.37
22	24.36
24	24.35
26	24.34
28	24.34
30	24.33



Test Setup Photo (Validation)

