


Note:
No CE mark for M8102AM



USA/Canada/EU/
Australia/New Zealand

Japan

Korea

Taiwan

Argentina

~~ITS #J45~~
(USA only)

WLAN #J35

SRR #J46

ITS #J47

電気式アナロイド型血圧計
IntelliVue MP2/X2
計量範囲 10~270mmHg
目 量 1mmHg
定格電圧 直流 10.8V
型式承認 第Q0818号

JAPAN ONLY

1709	Changed WLAN #J35 FCC ID and IC no.	01.03.17	Markus Rath
1641	Added symbol for Australia and label for Argentina	11.10.16	Markus Rath
1544	Changed Korea laser labeling WLAN & SRR, changed notified body from CE0366 to CE0123	27.10.15	Markus Rath
1540	Changed Russia symbol from GOST to EAC	06.10.15	Markus Rath
1531	Added Note	27.07.15	Markus Rath
Zchg.-Status DWG.- REV.	Änderungsbeschreibung / REVISION DESCRIPTION	Datum / DATE	Bestätigt / APPROVED
Maßneßteil Millimeter, falls nicht anders angegeben. / UNLESS OTHERWISE SPECIFIED DIMENSIONS ARE IN MILLIMETRES		Maßstab / SCALE	
Projektionsmethode / THIRD ANGLE PROJECTION 			
Nicht in der ausgedruckten Zeichnung messen. / DO NOT SCALE THIS DRAWING.		Datum/DATE 06.10.15	Name/NAME Markus Rath
- STRENG VERTRAULICH - WEITERGABE AN DRITTE UND VERVIELFÄLTIGUNGEN NUR MIT GENEHMIGUNG DES UHRHEBERS * SCHUTZVERMERK NACH ISO 16016 BEACHTEN *		Bearb./DRAWN Gepr./CHECKED	Benennung / TITLE Laser Labeling Rear MMS plus
PHILIPS MEDICAL SYSTEMS PROPRIETARY. This document contains confidential proprietary information that is PHILIPS PROPRIETARY. Do not disclose to or duplicate for others except as authorized by Philips Medical Systems.		Zeichnungsnummer / DRAWING NUMBER C-M3002-60223-1	
PHILIPS		Blatt / SHEET 1/1	
		Teilentr./PART NO.: M3002-60223	