



ortivus



MobiMedTM 300

Model:c



S/N. 0001240080

**Internally powered equipment.
Use only Ortivus power supply PS-1**

This device complies with part 15 of the FCC rules.
Operation is subject to the following two conditions:
(1) This device may not cause any harmful interference.
(2) This device must accept any interference received,
including interferece that may cause undesired operation.

FCC ID:PVJ-410000409



Conforms to
UL2601-1
Certified to
CAN/CSA
STD C22.2
No. 601.1-M90

Ortivus AB
Enhagsslingan 5
SE-187 40 Täby
SWEDEN



0413

IPX4



MADE IN SWEDEN