

REV	DESCRIPTION	D.C. NO.	DATE	DWN. BY	APPR. BY
A	PRODUCTION RELEASE	CJ446	06/29/01	HR	LH
B	SEE DCN	CP070	02/06/03	HR	KD
C	SEE DCN	C0021	4-30-03	HR	JD




PRECISIONTM
IMPLANTABLE PULSE GENERATOR KIT

REF



SC-1100



**MANUFACTURED BY
ADVANCED BIONICS® CORP.
12740 San Fernando Rd.
Sylmar, CA 91342 USA
(+1)818-362-7588**

ADVANCED BIONICS® SARL
76 Rue de Battenheim
68170 Rixheim, France
+ 33 (0) 3 89 65 98 00

SN



XXXXXXXXXX

STERILE LOT xxxxxxxx

FCC ID: Q4DSC-1100

6025421-001 Rev C

PRELIMINARY

1. USE LABEL #5063108-001
2. COLOR OF ALL LETTERS TO BE BLACK ON WHITE
3. INFORMATION IS ON COMPUTER IN LABELTRONIX EASYLABEL PROGRAM

NOTES: UNLESS OTHERWISE SPECIFIED:

DASH NO.	NEXT ASSY	MODEL NO.
APPLICATION		

UNLESS OTHERWISE SPECIFIED ALL DIMS ARE IN INCHES ALL DIMS APPLY AFTER FINISH		DRAIN A. CHEN		DATE 04/20/01		 ADVANCED BIONICS®		12740 San Fernando Road Sylmar, California 91342 1 800 678 2575	
TOLERANCES DECIMAL .XX ± .010 .XXX ± .005 .XXXX ± .0005		CHECKED K. CHINN		DATE 06/29/01		TITLE LABEL			
ANGLES ± 1° FRACTIONS ± 1/32		APPROVED L. HALBERG		DATE 06/29/01		SCS IPG TRAY			
SURFACE FINISH PER ANSI-B46-1 32/√		DO NOT SCALE DRAWING				SIZE B		DRAWING NUMBER 6025421	
BREAK ALL EDGES .005 MAX FILLET RADI .005 MAX		SCALE N/A		SHEET 1 of 1		REV C			

The label goes on the outside of the sterile tray. Since it's an implanted component, we don't have labels on the actual device.