


Blood pressure monitor  
part name: label,rating  
drwg: HEM-7371T1-AZ



\*Due to the size of product, the statement required by FCC 15.19(a) will be printed on users manual.

The FCC Identifier "FCC ID" is 4.0 point or larger.


The labels are laser-engraved, so the letters do not fade easily and the labels themselves do not peel off.

**BLOOD PRESSURE MONITOR**  
**MODEL: BP7460 [REF] HEM-7371T1-AZ**  
**AC ADAPTER: HEM-ADPTW5**  
**FCC ID: XXXXXXXXXX**  
 **OMRON HEALTHCARE Co., Ltd.**  
**KYOTO, 617-0002 JAPAN**  
**DISTRIBUTED BY: OMRON HEALTHCARE, INC.**  
**HOFFMAN ESTATES, IL 60169**  
**1-800-634-4350**  
**MADE IN VIETNAM**



FCC ID: XXX-XXXXXXX will be replaced by FCC ID:Q6ZHEM7382T1.

Blood pressure monitor  
part name: label,rating  
drwg: HEM-7372T1-AZAZ



\*Due to the size of product, the statement required by FCC 15.19(a) will be printed on users manual.

The FCC Identifier "FCC ID" is 4.0 point or larger.

The labels are laser-engraved, so the letters do not fade easily and the labels themselves do not peel off.


**BLOOD PRESSURE MONITOR/ TENSIO MÈTRE**  
**MODEL/ MODÈLE: BP5460 REF HEM-7372T1-AZAZ**  
**RATING/ CLASSMENT: DC6V 4.0W**  
**AC ADAPTER/ ADAPTEUR.A. : HEM-ADPTW5**  
**IC: XXXXX-XXXXXXXX**  
**FCC ID: XXXXXXXXXX**  
**CAN ICES-3(B)/ NMB-3(B)**  
 **OMRON HEALTHCARE Co., Ltd.**  
**KYOTO, 617-0002 JAPAN**  
**DISTRIBUTED BY: OMRON HEALTHCARE, INC.**  
**HOFFMAN ESTATES, IL 60169**  
**1-800-634-4350**  
**MADE IN VIETNAM**



IC:XXXXX-XXXXXXXX will be replaced by IC:10623A-HEM7382T1.

FCC ID: XXX-XXXXXXXX will be replaced by FCC ID:Q6ZHEM7382T1.

Blood pressure monitor  
part name: label,rating  
drwg: HEM-7381T1-AZ



\*Due to the size of product, the statement required by FCC 15.19(a) will be printed on users manual.


The FCC Identifier "FCC ID" is 4.0 point or larger.

The labels are laser-engraved, so the letters do not fade easily and the labels themselves do not peel off.



FCC ID: XXX-XXXXXXX will be replaced by FCC ID: Q6ZHEM7382T1.

Blood pressure monitor  
part name: label,rating  
drwg: HEM-7382T1-AZAZ



\*Due to the size of product, the statement required by FCC 15.19(a) will be printed on users manual.

The FCC Identifier "FCC ID" is 4.0 point or larger.

The labels are laser-engraved, so the letters do not fade easily and the labels themselves do not peel off.


**BLOOD PRESSURE MONITOR/ TENSIO MÈTRE**  
**MODEL/ MODÈLE: BP5465 REF HEM-7382T1-AZAZ**  
**RATING/ CLASSMENT: DC6V 4.0W**  
**AC ADAPTER/ ADAPTEUR.A. : HEM-ADPTW5**  
**IC: XXXXX-XXXXXXXX**  
**FCC ID: XXXXXXXXXX**  
**CAN ICES-3(B)/ NMB-3(B)**  
 **OMRON HEALTHCARE Co., Ltd.**  
**KYOTO, 617-0002 JAPAN**  
**DISTRIBUTED BY: OMRON HEALTHCARE, INC.**  
**HOFFMAN ESTATES, IL 60169**  
**1-800-634-4350**  
**MADE IN VIETNAM**



IC:XXXXX-XXXXXXXX will be replaced by IC:10623A-HEM7382T1.

FCC ID: XXX-XXXXXXXX will be replaced by FCC ID: Q6ZHEM7382T1.


Blood pressure monitor  
part name: label,rating  
drwg: HEM-7376T1-Z



\*Due to the size of product, the statement required by FCC 15.19(a) will be printed on users manual.


The FCC Identifier "FCC ID" is 4.0 point or larger.

The labels are laser-engraved, so the letters do not fade easily and the labels themselves do not peel off.



FCC ID: XXX-XXXXXXX will be replaced by FCC ID: Q6ZHEM7382T1.

Blood pressure monitor  
part name: label,rating  
drwg: HEM-7377T1-ZAZ



\*Due to the size of product, the statement required by FCC 15.19(a) will be printed on users manual.

The FCC Identifier "FCC ID" is 4.0 point or larger.

The labels are laser-engraved, so the letters do not fade easily and the labels themselves do not peel off.

**BLOOD PRESSURE MONITOR/ TENSIO MÈTRE**  
**MODEL/ MODÈLE: BP5360 REF HEM-7377T1-ZAZ**  
**RATING/ CLASSMENT: DC6V 4.0W**  
**AC ADAPTER/ ADAPTEUR.A. : HEM-ADPTW5**  
**IC: XXXXX-XXXXXXXX**  
**FCC ID: XXXXXXXXXX**  
**CAN ICES-3(B)/ NMB-3(B)**  
**OMRON HEALTHCARE Co., Ltd.**  
**KYOTO, 617-0002 JAPAN**  
**DISTRIBUTED BY: OMRON HEALTHCARE, INC.**  
**HOFFMAN ESTATES, IL 60169**  
**1-800-634-4350**  
**MADE IN VIETNAM**



IC:XXXXX-XXXXXXXXX will be replaced by IC:10623A-HEM7382T1.

FCC ID: XXX-XXXXXXXXX will be replaced by FCC ID: Q6ZHEM7382T1.