

Appendix 2 Photodocumentation

Description: Test setup - general view

Medizintechnik/Nachrichtentechnik/EMV

Department: EMC

Test report reference:
M/EMV-02/225

Page: 5 of 8

Date: 14. 10. 2002

checked by: 



Appendix 2 Photodocumentation

Description: Test setup for measurements below
30 MHz

Medizintechnik/Nach-
richtentechnik/EMV

Department: EMC

Test report reference:
M/EMV-02/225

Page: 6 of 8

Date: 14. 10. 2002

checked by: 



Appendix 2 Photodocumentation

Description: Test setup for measurements above
30 MHz

Medizintechnik/Nach-
richtentechnik/EMV


Department: EMC

Test report reference:
M/EMV-02/225

Page: 7 of 8

Date: 14. 10. 2002

checked by: 



Appendix 2 Photodocumentation

Description: Test setup with transponder card

Medizintechnik/Nachrichtentechnik/EMV

Department: EMC

Test report reference:
M/EMV-02/225

Page: 8 of 8

Date: 14. 10. 2002

checked by: *M*

