

Intica 7 HF-T QP

DF4/IS4



IS-1
IS4-LLLL
DF4-LLHH

REF 404630
Cat.No.

SN 12345678
Serial No.

PID 89
Prod.ID

Use By 2017-10-31

41 °F 113 °F
Ster.Temp.



Implantable Cardioverter-Defibrillator

Intica 7 HF-T QP

ProMRI®

PID 89
Prod.ID

Caution: Federal law (U.S.A.) restricts this device to sale by, or on the order of, a physician (or properly licenced practitioner)!


402-405 MHz
TP2



RA: IS-1 BI
LV: IS4-LLLL
RV: DF4-LLHH

FCC ID: QRTACHNT2

426020



For single use only, do not resterilize. Refer to technical manual for instructions.
ICD is shipped with all therapy disabled.

Manufacturer:
BIOTRONIK SE & Co. KG
Woermannkehre 1
12359 Berlin - Germany
Tel. +49 (0) 30 68905-0
Fax +49 (0) 30 6852804
sales@biotronik.com


Distributor:
BIOTRONIK, Inc.
6024 Jean Road
Lake Oswego, OR 97035-5369
(800) 547-0394 (24-hour)
(503) 635-9936 (Fax)



Intica 7 HF-T QP DF4/IS4


REF 404630
Cat.No.
SN 12345678
Serial No.
PID 89
Prod.ID

Use By 2017-10-31



Intica 7 HF-T QP

IS-1
IS4-LLLL
DF4-LLHH




DF4/IS4
ProMRI®

Implantable Cardioverter Defibrillator

REF 404630
Cat.No.
SN 12345678
Serial No.
PID 89
Prod.ID

STERILEEO
Sterilized with ethylene oxide



For single use only, do not resterilize.
Do not use if packaging is damaged.
Refer to technical manual for instructions.
ICD is shipped with all therapy disabled.

BIOTRONIK SE & Co. KG
Woermannkehre 1
12359 Berlin - Germany



(01)0403547914217917171031(2)12345678

SN 12345678
PID 89

SN 12345678
PID 89



c5034b2

Intica 7 HF-T QP DF4/IS4
SN 12345678
PID 89

Intica 7 HF-T QP DF4/IS4
SN 12345678
PID 89

Intica 7 HF-T QP DF4/IS4
SN 12345678
PID 89

Intica 7 HF-T QP DF4/IS4
SN 12345678
PID 89

Intica 7 HF-T QP DF4/IS4
REF 404630
SN 12345678
PID 89

Intica 7 HF-T QP DF4/IS4
REF 404630
SN 12345678
PID 89

Intica 7 HF-T QP DF4/IS4
REF 404630
SN 12345678
PID 89

Intica 7 HF-T QP DF4/IS4
REF 404630
SN 12345678
PID 89

Intica 7 HF-T QP DF4/IS4
REF 404630
SN 12345678
PID 89

Intica 7 HF-T QP

DF4/IS4 ProMRI® Data Sheet

IS-1
IS4-LLLL
DF4-LLHH

Implantable Cardioverter Defibrillator

Parameter	Values at Factory Release	Programmable Parameter Ranges
Bradycardia Therapy Mode	OFF	AAI(R); VOO; VVI(R); VVI-CLS; VDI(R); VDD(R); DOO; DDI(R); DDD(R); DDD-ADI(R); DDD-CLS; OFF
Basic Rate; min-1		30...[5]...100...[10]...160
Pulse Amplitude; V		0.5...[0.25]...4...[0.5]...6; 7.5
Pulse Width; ms		0.4; 0.5...[0.25]...1.5
A ATPs	OFF	OFF; BURST; RAMP; HF BURST
V ATPs	OFF	OFF; BURST; RAMP
Shock	OFF	OFF; DEFIBRILLATION with/without Confirmation
X-Ray Identification	NK	-

For single use only. Refer to Technical Manual for Instructions.

Distributor:
BIOTRONIK, Inc.
6024 Jean Road
Lake Oswego, OR 97035-5369
(800) 547-0394 (24-hour)
(503) 635-9936 (Fax)

BIOTRONIK SE & Co. KG
Woermannkehre 1
12359 Berlin - Germany
Tel. +49 (0) 30 68905-0
Fax +49 (0) 30 6852804
sales@biotronik.com



426020

