

ATTACHMENT P – SAR TEST SETUP PHOTOGRAPHS

■ Test Setup Photo



■ Test Setup Photo



■ Test Setup Photo



■ Test Setup Photo



■ Test Setup Photo




■ Test Setup Photo



■ Test Setup Photo





■ Test Setup Photo





■ Test Setup Photo

[835MHz/ 1900MHz Head Tissue Liquid Depth : 15cm]



■ Test Setup Photo

[835MHz/ 1900MHz Head Tissue Liquid Depth : 15cm]



■ Test Setup Photo

