


ATTACHMENT P – SAR TEST SETUP PHOTOGRAPHS


■ Test Setup Photo



■ Test Setup Photo



■ Test Setup Photo




■ Test Setup Photo




■ Test Setup Photo



■ Test Setup Photo




■ Test Setup Photo




■ Test Setup Photo




■ Test Setup Photo



■ Test Setup Photo




■ Test Setup Photo




■ Test Setup Photo



■ Test Setup Photo



■ Test Setup Photo



■ Test Setup Photo

[1900MHz Head Tissue Liquid Depth : 15cm]

