




EMC Laboratory
No.53-11, Tin-fu Tsun, Lin-kou Hsiang, Taipei
County, Taiwan R.O.C. Post Code:24443
Tel:+886-2-26092133 Fax:+886-2-26099303
Email:ttenc@ttenc.com.tw
Web:www.ttenc.com.tw




EMC Laboratory
No.53-11, Tin-fu Tsun, Lin-kou Hsiang, Taipei
County, Taiwan R.O.C. Post Code:24443
Tel:+886-2-26092133 Fax:+886-2-26099303
Email:ttenc@ttenc.com.tw
Web:www.ttenc.com.tw

Data#: 5 File#: C:\K-jump-2004-03TX.emi



Site : Anechoic Chamber
Condition : FCC CLASS-B 3m BBA9106(A3L) HORIZONTAL
EUT : Clinic Patch Thermometer
: (Stand)
MEMO : M/N:KD-2100
POWER : DC 6V
ENVIRONMENT : 17°C/46%

Data#: 6 File#: C:\K-jump-2004-03TX.emi



Site : Anechoic Chamber
Condition : FCC CLASS-B 3m BBA9106(A3L) VERTICAL
EUT : Clinic Patch Thermometer
: (Stand)
MEMO : M/N:KD-2100
POWER : DC 6V
ENVIRONMENT : 17°C/46%




EMC Laboratory
No.53-11, Tin-fu Tsun, Lin-kou Hsiang, Taipei
County, Taiwan R.O.C. Post Code:24443
Tel:+886-2-26092133 Fax:+886-2-26099303
Email:ttenc@ttenc.com.tw
Web:www.ttenc.com.tw




EMC Laboratory
No.53-11, Tin-fu Tsun, Lin-kou Hsiang, Taipei
County, Taiwan R.O.C. Post Code:24443
Tel:+886-2-26092133 Fax:+886-2-26099303
Email:ttenc@ttenc.com.tw
Web:www.ttenc.com.tw

Data#: 12 File#: C:\K-jump-2004-03TX.emi



Site : Anechoic Chamber
Condition : FCC CLASS-B 3m UHALP 9108-A 0138 HORIZONTAL
EUT : Clinic Patch Thermometer
: (Stand)
MEMO : M/N:KD-2100
POWER : DC 6V
ENVIRONMENT : 17°C/46%

Data#: 11 File#: C:\K-jump-2004-03TX.emi



Site : Anechoic Chamber
Condition : FCC CLASS-B 3m UHALP 9108-A 0138 VERTICAL
EUT : Clinic Patch Thermometer
: (Stand)
MEMO : M/N:KD-2100
POWER : DC 6V
ENVIRONMENT : 17°C/46%




EMC Laboratory
No.53-11, Tin-fu Tsun, Lin-kou Hsiang, Taipei
County, Taiwan R.O.C. Post Code:24443
Tel:+886-2-26092133 Fax:+886-2-26099303
Email:ttenc@ttenc.com.tw
Web:www.ttenc.com.tw




EMC Laboratory
No.53-11, Tin-fu Tsun, Lin-kou Hsiang, Taipei
County, Taiwan R.O.C. Post Code:24443
Tel:+886-2-26092133 Fax:+886-2-26099303
Email:ttenc@ttenc.com.tw
Web:www.ttenc.com.tw

Data#: 4 File#: C:\K-jump-2004-03TX.emi



Site : Anechoic Chamber
Condition : FCC CLASS-B 3m BBA9106(A3L) HORIZONTAL
EUT : Clinic Patch Thermometer
: (Side)
MEMO : M/N:KD-2100
POWER : DC 6V
ENVIRONMENT : 17°C/46%

Data#: 3 File#: C:\K-jump-2004-03TX.emi



Site : Anechoic Chamber
Condition : FCC CLASS-B 3m BBA9106(A3L) VERTICAL
EUT : Clinic Patch Thermometer
: (Side)
MEMO : M/N:KD-2100
POWER : DC 6V
ENVIRONMENT : 17°C/46%




EMC Laboratory
No.53-11, Tin-fu Tsun, Lin-kou Hsiang, Taipei
County, Taiwan R.O.C. Post Code:24443
Tel:+886-2-26092133 Fax:+886-2-26099303
Email:ttenc@ttenc.com.tw
Web:www.ttenc.com.tw




EMC Laboratory
No.53-11, Tin-fu Tsun, Lin-kou Hsiang, Taipei
County, Taiwan R.O.C. Post Code:24443
Tel:+886-2-26092133 Fax:+886-2-26099303
Email:ttenc@ttenc.com.tw
Web:www.ttenc.com.tw

Data#: 9 File#: C:\K-jump-2004-03TX.emi



Site : Anechoic Chamber
Condition : FCC CLASS-B 3m UHALP 9108-A 0138 HORIZONTAL
EUT : Clinic Patch Thermometer
: (Side)
MEMO : M/N:KD-2100
POWER : DC 6V
ENVIRONMENT : 17°C/46%

Data#: 10 File#: C:\K-jump-2004-03TX.emi



Site : Anechoic Chamber
Condition : FCC CLASS-B 3m UHALP 9108-A 0138 VERTICAL
EUT : Clinic Patch Thermometer
: (Side)
MEMO : M/N:KD-2100
POWER : DC 6V
ENVIRONMENT : 17°C/46%




EMC Laboratory
No.53-11, Tin-fu Tsun, Lin-kou Hsiang, Taipei
County, Taiwan R.O.C. Post Code:24443
Tel:+886-2-26092133 Fax:+886-2-26099303
Email:ttenc@ttenc.com.tw
Web:www.ttenc.com.tw




EMC Laboratory
No.53-11, Tin-fu Tsun, Lin-kou Hsiang, Taipei
County, Taiwan R.O.C. Post Code:24443
Tel:+886-2-26092133 Fax:+886-2-26099303
Email:ttenc@ttenc.com.tw
Web:www.ttenc.com.tw

Data#: 1 File#: C:\K-jump-2004-03TX.emi



Site : Anechoic Chamber
Condition : FCC CLASS-B 3m BBA9106(A3L) HORIZONTAL
EUT : Clinic Patch Thermometer
(Lying)
MEMO : M/N:KD-2100
POWER : DC 6V
ENVIRONMENT : 17°C/46%

Data#: 2 File#: C:\K-jump-2004-03TX.emi



Site : Anechoic Chamber
Condition : FCC CLASS-B 3m BBA9106(A3L) VERTICAL
EUT : Clinic Patch Thermometer
(Lying)
MEMO : M/N:KD-2100
POWER : DC 6V
ENVIRONMENT : 17°C/46%




EMC Laboratory
No.53-11, Tin-fu Tsun, Lin-kou Hsiang, Taipei
County, Taiwan R.O.C. Post Code:24443
Tel:+886-2-26092133 Fax:+886-2-26099303
Email:ttenc@ttenc.com.tw
Web:www.ttenc.com.tw




EMC Laboratory
No.53-11, Tin-fu Tsun, Lin-kou Hsiang, Taipei
County, Taiwan R.O.C. Post Code:24443
Tel:+886-2-26092133 Fax:+886-2-26099303
Email:ttenc@ttenc.com.tw
Web:www.ttenc.com.tw

Data#: 6 File#: C:\K-jump-2004-03TX.emi



Site : Anechoic Chamber
Condition : FCC CLASS-B 3m UHALP 9108-A 0138 HORIZONTAL
EUT : Clinic Patch Thermometer
(Lying)
MEMO : M/N:KD-2100
POWER : DC 6V
ENVIRONMENT : 17°C/46%

Data#: 7 File#: C:\K-jump-2004-03TX.emi



Site : Anechoic Chamber
Condition : FCC CLASS-B 3m UHALP 9108-A 0138 VERTICAL
EUT : Clinic Patch Thermometer
(Lying)
MEMO : M/N:KD-2100
POWER : DC 6V
ENVIRONMENT : 17°C/46%