



REVISIONS				
ZONE	LTR	DESCRIPTION	DATE	APPROVED
	x3	EO #		



WRITE IN 6 DIGIT NUMBER DISPLAYED
ON TEST RECEIVER USING FINE POINT
PERMANENT MARKER, THEN INSTALL
LABEL ON POUR SPOUT.

#	DESCRIPTION	CODE IDENT
2800	Vital Link SPRCV-2 Label FCC	1
75008	Vital Link SPRCV-2 Cork large	2
2400	S. A. PT-200 Pour Spout	3
3900	Procedure - PT-200 Serial Number ID Check and Label Install	

NOTE: TUBE IS INSTALLED AT VITAL LINK

UNLESS OTHERWISE SPECIFIED
DIMENSIONS ARE IN INCHES
TOLERANCES ARE:
FRACTIONS DECIMALS ANGLES
 $\pm 1/64$.XX $\pm .010$ $\pm$
.XXX $\pm .005$

RF Applications, Inc.

ASSY. PT-200 POUR SPOUT
W/LARGE CORK

NOV 14 2002
CHECK PRINT
APPROVED FOR RELEASE
DATE 11/14/02 approved all

DRAWN BY
CA
DATE
10-29-02

CHECKED BY
DATE

APPROVED BY
DATE

REVISION
DATE

SIZE
B
PROJECT NO.
02E011

MATERIAL
FINISH

SHRINKAGE
COLOR

PART NO.
33B37210

REV.
x3

DO NOT SCALE DRAWING

SCALE 1=1

SHEET 1 OF 1

9 2002