



REVISIONS				
ZONE	LTR	DESCRIPTION	DATE	APPROVED
	x3	EO #		



WRITE IN 6 DIGIT NUMBER DISPLAYED ON TEST RECEIVER USING FINE POINT PERMANENT MARKER, THEN INSTALL LABEL ON POUR SPOUT.

#	DESCRIPTION	CODE IDENT
2800	Vital Link SPRCV-2 Label FCC	1
75008	Vital Link SPRCV-2 Cork large	2
2400	S. A. PT-200 Pour Spout	3
3900	Procedure - PT-200 Serial Number ID Check and Label Install	

NOTE: TUBE IS INSTALLED AT VITAL LINK

UNLESS OTHERWISE SPECIFIED
DIMENSIONS ARE IN INCHES
TOLERANCES ARE:
FRACTIONS DECIMALS ANGLES
 $\pm 1/64$.XX $\pm .010$ $\pm$
.XXX $\pm .005$

RF Applications, Inc.

DRAWN BY
CA
DATE
10-29-02

ASSY. PT-200 POUR SPOUT
W/LARGE CORK

CHECKED BY
DATE

APPROVED BY
DATE

MATERIAL

FINISH

REVISION
DATE

SHRINKAGE

COLOR

SIZE
B

PROJECT NO.
02E011

PART NO.

33B37210

REV.
x3

DO NOT SCALE DRAWING

SCALE 1=1

SHEET 1 OF 1

NOV 14 2002
CHECK PRINT

APPROVED FOR RELEASE

DATE 11/14/02 approved *all*

9 2002

FCC ID: TGJ-SPRCV1-F4

IC: 6245A-SPRCV1F4

S/N



Product Code:

SPRCV-4

Capton, Inc

USA