

# ID Label Location & Sample

**LABEL LOCATION:** Label will be attached to the Rear of the product



**LABEL Drawing:**

|                                                                                                                                                                                                                                        |                                                                                     |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|
| BLOOD GLUCOSE MONITORING SYSTEM                                                                                                                                                                                                        |                                                                                     |
| • MODEL : IGM-1007A                                                                                                                                                                                                                    |                                                                                     |
| • RATING : DC 3V $\frac{1}{2}$ , 17mA                                                                                                                                                                                                  |                                                                                     |
| REF                                                                                                                                                                                                                                    |                                                                                     |
| SN                                                                                                                                                                                                                                     |                                                                                     |
|  <b>OSANG Healthcare Co., Ltd.</b><br>132, Anyangcheondong-ro, Dongan-gu,<br>Anyang-si, Gyeonggi-do, Korea<br>www.osanghc.com                       |                                                                                     |
|  <b>Obelis S.A</b><br>Bd. General Wahis 53, 1030 Brussels, Belgium<br>Tel : +(32) 2.732.59.54 / Fax : +(32) 2.732.60.03<br>E-mail : mail@obelis.net |                                                                                     |
| IVD                                                                                                                                                                                                                                    | SELF-TESTING                                                                        |
|                                                                                                                                                     |  |
|                                                                                                                                                     |  |
|  <b>0123</b>  R-C-IPM-IGM-1007A                                  |                                                                                     |
| FCC ID: WSX-IGM-1007A                                                                                                                                                                                                                  |                                                                                     |